




Philips Sonicare AirFloss  
Rechargeable powered  
interdental cleaner

w/3-series handle  
w/ 1 nozzle and 1 brush head



HX8274/20

## Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard-to-reach areas.

### Improves oral health

- Helps improve gum health in just two weeks
- Helps prevent cavities between teeth
- Remove up to 7x more plaque with our InterCare brush head\*

### Achieve a 30-second clean

- Cleans your entire mouth in just 30 seconds
- Point. Press. Clean!

### Easy to use

- Air and micro-droplet technology
- Guidance tip ensures correct placement

### Helps build a healthy habit

- An easy way to improve cleaning between teeth

**PHILIPS**  
**sonicare**


# Highlights

## Air and micro-droplet tech



AirFloss utilises unique air and micro-droplet technology to send a burst of pressurised air and water or mouthwash between your teeth. The one-touch activation button is as simple as a mouse click.

## An easy way to start




Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

## 30-Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

## Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks.

## Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

## Guidance tip



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

## Helps prevent cavities



By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

## Remove up to 7 x more plaque



Click on our InterCare brush head to improve gum health in just 2 weeks. Extra-long bristles help remove more plaque in hard-to-reach spots and between teeth for healthy gums.

# Specifications

## Power

- Voltage: Multi-voltage charger

## Technical specifications

- Operating time: Replace nozzle every 6 months
- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

## Design specifications

- Colour(s): White with green accents

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Handle: Slim ergonomic design

- Battery indicator: Flashing indicates low battery
- Battery Life: Up to 3 weeks when fully charged
- Nozzle attachment: Easily snaps on and off

## Items included

- Charger: 1 multi-voltage charger
- AirFloss Handle: 1
- Handle: 1 3-Series
- AirFloss Nozzle: 1
- Brush head: 1 InterCare standard

## Cleaning performance

- Auto-burst feature: Cleans entire mouth in 30 secs
- Health benefits: Improves gum health in only two weeks
- Nozzles: Best results change every 6m
- Performance: Removes up to 7x more plaque\*



Issue date 2023-05-19

Version: 5.0.2

EAN: 87 10103 82327 8

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* than manual tooth-brushing alone