

PHILIPS

sonicare

Sonic electric
toothbrush

DiamondClean

5 modes

2 brush heads

Glass charger

USB travel charger

HX9332/12



Whiter, healthier teeth for life

Best Philips Sonicare whitening sonic toothbrush

Sonicare's best whitening in our most elegant Philips Sonicare electric toothbrush. Switch to Sonicare.

Proven to improve oral health

- Improves gum health in only two weeks
- Philips Sonicare toothbrush helps whiten teeth
- Removes up to 7 x more plaque than a manual toothbrush
- Whitens teeth more than 2x better than a manual toothbrush
- Sensitive Mode: Gentle teeth and gum cleaning

Designed around you

- Dual Charging System: charging glass and travel case.

Guides you to follow dental recommendations


- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

Provides a superior clean

- DiamondClean toothbrush head for Sonicare's best whitening
- Philips Sonicare's advanced sonic technology


Highlights

Removes up to 7 x more plaque



Removes up to 7 times more plaque from hard-to-reach places than a manual toothbrush.

Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

Whitens teeth twice as well



This Philips Sonicare toothbrush whitens teeth more than 2x better than a manual toothbrush in just 1 week. Removes up to 100% more stains for whiter teeth in just 1 week.

DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles to effectively yet


gently remove surface stains and scrape away plaque. Removes stain build up from coffee, tea, red wine, cigarettes and tobacco in just 1 week.

Dual Charging System



1) Charger glass - Simply place your Philips Sonicare electric toothbrush in the glass to charge. Can also be used to rinse your mouth after brushing. 2) Charging travel case - Connects to your laptop via USB or can be charged by wall socket. Allows you to take your DiamondClean anywhere and everywhere.

Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Smartimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Helps whiten teeth



This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

The only true sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.



Specifications

Power

Voltage: 110–220 V

Technical specifications

Operating time (full to empty): Up to 3 weeks**

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Colour: Ceramic white

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Illuminated icon indicates battery life

Display: Illuminated display

Brushing time: Up to 3 weeks**

Items included

Travel case: USB travel charger

Brush heads: 2 DiamondClean standard

Handles: 1 DiamondClean

Glass charger: 1

Cleaning performance

Health benefits: Improves gum health in only two weeks

Speed: Up to 62,000 brush movements/min

Performance: Removes up to 7x more plaque*

Whitening benefits: Whitens teeth twice as well*

Timer: Quadpacer and Smartimer

Modes

Clean: For exceptional everyday cleaning

Gum Care: Gently massages gums

White: Removes surface stains

Polish: Brighten and polish your teeth

Sensitive: Gentle teeth and gum cleaning

