

PHILIPS

sonicare

Sonic electric
toothbrush

FlexCare

3 modes

2 brush heads

UV brush head sanitiser

HX6932/10



Advanced Cleaning

Toothbrush for gentle cleaning of teeth and gums

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare electric toothbrush that adapts to your oral care needs.

Proven to improve oral health

- Improves gum health in only two weeks
- Philips Sonicare toothbrush helps whiten teeth
- Sensitive Mode: Gentle teeth and gum cleaning

Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

Provides a superior clean

- Sonicare dynamic sonic action drives fluid between teeth
- Removes up to 2x more plaque than a manual toothbrush

Offers a customised brushing experience

- Easy-start programme to acclimate to the Sonicare experience
- Three different modes for a better cleaning experience


Highlights

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

Smartimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Quadpacer




30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience

Three different cleaning modes




The three modes include a Clean mode for an outstanding clean in 2 minutes. A Sensitive mode for gentle yet thorough cleaning of the gums and a Massage mode to stimulate the gums.

Helps whiten teeth



This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums



Specifications

Design and finishing

Colour: White and green

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: Lithium ION

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design, Rubber grip for easy handling

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

Items included

Travel case: 1 soft

Brush heads: 1 ProResults standard, 1 ProResults compact

Charger: 1

Handles: 1 FlexCare

UV sanitiser

Cleaning performance

Health benefits: Helps improve gum health

Speed: Up to 62,000 brush movements/min

Performance: Removes up to 83% of plaque in hard-to-reach areas

Whitening benefits: Whitens teeth by up to 2 shades

Timer: Quadpacer and Smartimer

Modes

Clean: For exceptional everyday cleaning

Massage: Invigorating massage

Sensitive: Gentle teeth and gum cleaning

2 Routines: Go Care and Max Care

