

PHILIPS

sonicare

Sonic electric
toothbrush

ProtectiveClean
4100

Smart brush head sensors

1 cleaning mode

1 x BrushSync feature



HX6888/18



Cleaner teeth. Gently does it.

Remove up to 7x more plaque than a manual

Feel the difference of a gentle clean with our Pressure Sensor while removing up to 7x more plaque than a manual toothbrush

Provides a superior clean

- Remove up to 7x more plaque with our InterCare brush head*

Innovative technology

- A gentle nudge to go nice and easy
- Patients gets the most from their brush heads
- Connects smart brush handle and smart brush heads
- Philips Sonicare's advanced sonic technology

Improves patients' oral health

- Fillings, crowns, or sensitive areas? Not to worry.

Designed around your patient

- A great way to get used to your new toothbrush
- Encouragement to brush thoroughly
- Just the ticket for patients on the go


Highlights

Remove up to 7x more plaque



Click on our InterCare brush head to improve gum health in just 2 weeks. Extra-long bristles help remove more plaque in hard-to-reach spots and between teeth for healthy gums.

Too much pressure?



Patients may not notice if they are brushing too hard, but their toothbrush will. If patients start taking it out on their gums, the toothbrush will make a pulsing sound as a reminder to ease off on the pressure.

Easy-start program



Starting something new can take getting used to. That's why our easy-start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

Safe and gentle



You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns, and veneers, and helps prevent cavities and improve gum health.

BrushSync technology



A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

Handy travel case and charger




Our premium travel case lets patients store their toothbrush hygienically, while our compact travel charger keeps them topped up while on the go. Patients will get two weeks of regular use from a single full charge, but the charger's just the ticket for longer trips.

Philips Sonicare technology




Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Smart replacement reminders



All brush heads wear out over time. But we've got that covered too: our smart brush-recognition technology tracks how long patients have been using their brush head for, and how hard they've been brushing. When it's time to replace it, a light on the handle lets them know. That way, they can be sure that their brush head is consistently doing a stellar job.

Helpful timers



Our QuadPacer lets patients know when they've spent just the right amount of time cleaning each part of their mouth, while our Smartimer tells them when they've brushed for the recommended two minutes.

Specifications

Design and finishing

Color: White and Mid Blue

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: Lithium ION

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

Handle compatibility: Easy click-on brush

heads

Brushing time: Up to 2 weeks

Replacement reminder: To always ensure best results, reminder icon lights up

Items included

Charger: 1

Brush heads: 2 i InterCare

Travel case: 1

Handle: 1 ProtectiveClean

Cleaning performance

Health benefits: Improves gum health in only two weeks, Helps reduce cavities

Performance: Removes up to 7x more plaque*

Speed: Up to 62000 brush movement/min

Timer: Quadpacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday clean

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

