

Philips Sonicare
TongueCare+
TongueCare+ spray kit

1 tongue brush
1 tongue spray
BrushSync mode pairing

HX8071/33



Clean tongue, fresh breath

A Philips Sonicare clean, for your tongue

Click. Spray. Clean. Our breakthrough smart TongueCare+ tongue brush and spray instantly transform your Philips Sonicare toothbrush into a powerful fresh breath solution.

Your tongue gets a Sonicare clean. You get long-lasting fresh breath.

An easy addition to your oral care routine

- From toothbrush to tongue cleaner in just one click
- Kill bad-breath bacteria and neutralise odours instantly
- Easily reach every part of your tongue

Innovative technology

- Automatically pairs with your toothbrush handle
- 240 MicroBristles for a gentle deep clean
- A Philips Sonicare clean for your tongue

PHILIPS
sonicare

TongueCare+ spray kit
1 tongue brush 1 tongue spray, BrushSync mode pairing

HX8071/33

Highlights

BreathRx tongue spray




Smelly VSC gases are no match for our BreathRx tongue spray, which is proven to kill odour-causing bacteria and neutralise bad breath instantly. The special formula coats your tongue and keeps it lubricated for a gentle cleaning experience with the TongueCare+ tongue brush. The cool mint flavour leaves your mouth feeling fresh and clean.

Easy to click on




The TongueCare+ tongue brush quickly transforms your Philips Sonicare toothbrush into a sonic-powered tongue cleaner. Each tongue brush clicks onto your Philips Sonicare toothbrush handle, just like a standard brush head. Simple to use, replace and clean, these innovative tongue brushes are compatible with all Philips Sonicare click-on toothbrush handles, so it couldn't be easier to add tongue cleaning to your daily oral care routine. It fits all Philips Sonicare toothbrush handles except PowerUp Battery and Essence.

Get your tongue truly clean



To get your tongue truly clean, you need to be able to reach everywhere. Thanks to its ergonomic shape, the TongueCare+ tongue brush lets you do this in comfort. The small and flexible head makes it easy to clean away bad-breath bacteria in every nook and cranny.

Smart mode-pairing



Our smart TongueCare+ tongue brush has a microchip inside that communicates with the toothbrush handle. So when you click it on, the handle automatically recognises it and selects the appropriate mode and intensity level to optimise your tongue-cleaning performance. All you have to do is to press the power button.

Microbristles



With 240 rubber MicroBristles, the TongueCare+ tongue brush is specifically designed to clean the soft and porous surface of your tongue. The flexible MicroBristles contour to and gently clean around the ridges and grooves, to remove bacteria build-up and drive our bacteria-killing Philips Sonicare BreathRx tongue spray deeper. The compact shape allows you to comfortably clean the entire tongue.

Sonic technology



TongueCare+ is the only tongue cleaner with Philips Sonicare strength, to clean away bad-breath bacteria with every stroke. 31,000 vibrations per minute break up tongue debris, sweep away bad-breath bacteria and drive the bacteria-killing ingredients of the Philips Sonicare BreathRx tongue spray deeper. Your tongue will feel cleaner and your breath will stay fresh longer.

Specifications

Design and finishing

- Colour: Black

Material

- Elastomer: Latex-free

Compatibility

- Brush head system: Click-on
- BrushSync mode pairing
- Suitable for these models: 2 Series plaque defence, 2 Series plaque defence, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean, Essence+, FlexCare, FlexCare Platinum, FlexCare Platinum Connected, FlexCare+, HealthyWhite, HealthyWhite+, PowerUp

Items included

- Brush heads: 1 TongueCare+ tongue brush
- Antibacterial tongue spray: 2 oz

Ingredients

- Tongue Spray Ingredients: Aqua, Propylene Glycerol, Aroma, Zinc Gluconate, Cetylpyridinium Chloride, Eucalyptus, Citriodora Oil, Limonene, CI42090, PEG-12 Dimethicone, PEG-40 Hydrogenated Castor Oil, Thymol, Sodium Saccharin

Use by

- After opened: 6 months
- Tongue spray shelf life: 36 months



Issue date 2022-06-04

Version: 3.0.2

EAN: 08 71010 38055 40

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com