

PHILIPS

sonicare

Sonic electric
toothbrush

5 Series gum health

2 modes

6 customizable settings

1 brush head

HX8931/10



Healthier gums in just 2 weeks

Reduces redness, swelling and bleeding gums

Provides plaque removal and gum care to help prevent and reduce the early signs of gum disease while offering your most comfortable brushing experience. Clinically proven to improve gum health 3 ways: reduces redness, swelling, bleeding.

Proven to improve oral health

- Proven to improve gum health 3 ways
- Better than a manual toothbrush for gingivitis
- Removes up to 6X more plaque along the gumline

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Offers a customized brushing experience

- Clean and Gum Plus mode
- 6 customizable settings

Provides a superior and gentle clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Safe on orthodontics, dental work and veneers

Highlights

Improves gum health 3 ways



Clinically proven to improve gum health three ways: reduces redness, swelling, bleeding (sites)

Prevents & reduces gingivitis



With superior plaque removal along the gum line, this toothbrush can help reduce gingivitis up to 100% better than a manual toothbrush in just two weeks

Safe and gentle



The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worry.

Up to 6x better plaque removal



Removes up to 6X more plaque along the gumline than a manual toothbrush.

Sonic technology




Our sonic technology creates a unique dynamic cleaning action that drives fluid deep between teeth and along the gumline for a thorough-yet-gentle clean

6 customizable settings



Philips Sonicare toothbrush allows you to mix and match modes (Clean mode, Gum Plus mode) and 3 intensities for maximum comfort and your best clean.

Clean and Gum Plus mode




2 minute Clean mode offers the ultimate in plaque removal (Dental professional recommended brushing time). Gum Plus mode stimulates and gently massages gums.

SmartTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth



Specifications

Design and finishing

Color: Frost blue

Power

Voltage: 110-220 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 3 weeks**

Battery type: Lithium ION

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Items included

Handles: 1 5 Series gum health

Brush heads: 1 G2 Optimal Gum Care

Charger: 1

Cleaning performance

Health benefits: Improves gum health in only two weeks

Speed: Up to 31000 brush movement/min

Performance: Removes up to 6X more plaque*

Timer: Quadpacer and SmarTimer

Modes

3 intensities: Low, Medium, High

Clean: For exceptional everyday clean

Gum Plus: Stimulates gums

