

PHILIPS

sonicare

Sonic electric
toothbrush

FlexCare

3 modes

2 brush heads

HX6942/04



Redefine clean

with FlexCare features

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare that adapts to your oral care needs.

Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

ProResults brush head

- Designed for a better fit with your mouth

FlexCare features

- 3 flexible brushing modes, 2 personalised cleaning routines
- Max care: 3 minutes for an extended care routine
- Go Care: 1-minute brushing cycle for a quick clean

Highlights

ProResults brush head



ProResults brush head's contour fit, wider bristle range and broader sweeping motion provide better tooth coverage.

FlexCare features



3 flexible brushing modes, 2 personalised cleaning routines.

Go Care


One-minute brushing cycle in the Clean mode, for a quick clean.



Max care

3 minutes for an extended care routine.

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Smartimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Design and finishing

Colour: White and green

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: Lithium ION

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design, Rubber grip for easy handling

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

Items included

Brush heads: 1 ProResults standard, 1 ProResults compact

Charger: 1

Handles: 1 FlexCare

Cleaning performance

Health benefits: Helps improve gum health

Speed: Up to 62,000 brush movements/min

Performance: Removes up to 83% of plaque in hard-to-reach areas

Whitening benefits: Whitens teeth by up to 2 shades

Timer: Quadpacer and Smartimer

Modes

Clean: For exceptional everyday cleaning

Massage: Invigorating massage

Sensitive: Gentle teeth and gum cleaning

2 Routines: Go Care and Max Care

