

# PHILIPS

## sonicare

Sonic electric  
toothbrush

HealthyWhite

3 Modes

2 brush heads

HX6732/45



## Whiter, healthier teeth

### The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Sonicare HealthyWhite electric toothbrush is proven to remove everyday stains in just two weeks, with the regular use of the Clean and White Mode.

#### **Proven to improve oral health**

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology

#### **Guides you to follow dental recommendations**

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

#### **Provides a superior clean**

- Sonicare dynamic sonic action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

#### **Offers a customised brushing experience**

- Easy-start programme to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2-minute mode
- Sensitive Mode: Gentle teeth and gum cleaning

# Highlights

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

## Angled brush head neck



The uniquely angled neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Naturally whiter teeth



The dynamic sonic action of this toothbrush and the direct broader contact with each tooth, has been proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## Smartimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer




30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

## Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience

## Clean and White Mode




2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.\*

## Clean Mode

For exceptional everyday cleaning. Helps maintain whitening treatments.

## Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums



# Specifications

## Power

Voltage: 110–220 V

## Technical specifications

Operating time (full to empty): Up to 3 weeks\*\*

Battery: Rechargeable

Battery type: Lithium ION

## Design and finishing

Colour: Frost White

## Service

Warranty: 2-year limited warranty

## Ease of use

Handle: Slim ergonomic design, Rubber grip for easy handling

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 3 weeks\*\*

## Items included

Travel case: 1

Brush heads: 1 ProResults standard, 1 DiamondClean standard

Charger: 1

Handles: 1 HealthyWhite

## Cleaning performance

Health benefits: Helps improve gum health

Speed: Up to 62,000 brush movements/min

Performance: Removes up to 90% of stains\*

Whitening benefits: Whitens teeth by up to 2 shades

Timer: Quadpacer and Smartimer

## Modes

Clean: For exceptional everyday cleaning

Sensitive: Gentle teeth and gum cleaning

Clean and White: 2 min Clean + 30 sec White

