

PHILIPS

sonicare

Sonic electric
toothbrush

HealthyWhite

2 modes

1 brush head

HX6711/06



Whiter, healthier teeth

The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Philips Sonicare HealthyWhite electric toothbrush is proven to remove everyday stains in just two weeks, with the use of the Clean and White Mode.

Proven to improve oral health

- Philips Sonicare toothbrush helps whiten teeth

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Philips Sonicare's advanced sonic technology
- Removes up to 7 x more plaque than a manual toothbrush
- DiamondClean brush head offering the best whitening

Offers a customised brushing experience

- Easy-start programme to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains

Highlights

Helps whiten teeth



This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*


DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond-shaped, medium-stiffness bristles to

effectively yet gently scrape away plaque. This Philips Sonicare electric toothbrush provides exceptional cleaning and whiter teeth compared to a manual toothbrush.

Removes up to 7 x more plaque




Removes up to 100% more plaque from hard to reach places than a manual toothbrush.

The only true sonic technology




Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Smartimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time



Specifications

Power

Voltage: 110–220 V

Technical specifications

Operating time (full to empty): Up to 3 weeks**

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Colour: Glacier Green

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 3 weeks**

Items included

Brush heads: 1 DiamondClean standard

Charger: 1

Handles: 1 HealthyWhite

Cleaning performance

Health benefits: Helps improve gum health

Speed: Up to 62,000 brush movements/min

Performance: Removes up to 7x more plaque*

Whitening benefits: Helps naturally remove stains, Whitens teeth by up to 2 shades

Timer: Quadpacer and Smartimer

Modes

Clean: For exceptional everyday cleaning

Clean and White: 2 min Clean + 30 sec White

