

PHILIPS

AVENT

Breastfeeding
Support Kit

Natural

SCD292/00



Breastfeeding made easy

Pump, store, feed, care

The Natural Breastfeeding Support Kit is the ideal set for mums who have just started expressing milk. The set contains all the products you need to express, store and feed with your breast milk and to take care of your breasts.

Specifications

Country of origin

England

Development stages

Stage: Pregnancy, 0-6 months

Material

Bottle: Polypropylene, BPA-free*

Breast pump: BPA-free*, Polypropylene

Teat: Silicone, BPA-free*

Breast pads: Silky soft top-sheet, Natural breathable material, Dermatologically tested

Storage cups: Polypropylene, BPA-free*

Nipple protector: Silicone, BPA Free*

What is included

Milk storage cup lid: 10 pcs

Milk storage cup adapter: 2 pcs

Sealing disc for milk storage: 1 pcs

Extra soft Newborn Flow teat: 1 pcs

Extra soft Slow Flow teat: 1 pcs

Design

Bottle design: Ergonomic shape, Wide neck

Breast pump design: Compact design

What is included

Natural bottle 125 ml/4 oz: 1 pcs

Natural bottle 9 oz: 1 pcs

Breast pump body: 1 pcs

Standard-size cushion: 1 pcs

Base unit incl. tubing: 1 pcs

Milk storage cup (180 ml/6 oz): 10 pcs

Travel cover: 1 pcs

Guide to breastfeeding DVD: 1 pcs

Disposable day breast pads: 12 pcs

Disposable night breast pads: 6 pcs

Standard Nipple Protector: 2 pcs

Bra clip: 1 pcs

Moisturising Nipple Cream: 15 ml

Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Settings: 3 Expression settings, 1 Stimulation mode

