

PHILIPS

sonicare

Sonic electric
toothbrush

3100 series

Sonic technology

Pressure sensor

QuadPacer and SmarTimer

Slim ergonomic design



HX3671/54



Goodbye manual toothbrush. Hello Sonic technology.

Up to 3x more plaque removal*

Sonic technology combined with our brush action gently removes 3x more plaque than a manual toothbrush.

Designed around you

- 14-day long lasting battery life
- Ergonomic design makes the toothbrush easy to hold and use
- Always know when to replace your brush heads

Optimize your brushing

- Optimize your brushing with SmarTimer and QuadPacer

Easy transition to Sonicare

- Easy-Start to ease the transition

Protects your gums


- Built-in pressure sensor protects your teeth & gums

Provides superior clean

- Removes up 3x more plaque*vs a manual toothbrush
- Our unique technology brings you a powerful yet gentle clean

Highlights

Advanced Sonicare Technology



Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two month's worth of manual brushing in just 2 minutes.** 31000 brush strokes per minute gently clean your teeth, break up plaque, and sweep it away for an exceptional daily clean.

Removes up to 3X more plaque*



Sonicare electric toothbrush with advanced sonic technology is clinically proven to remove upto 3X more plaque from your teeth and along the gum line while protecting your gums

Protects your gums



Built-in pressure sensor automatically detects pressure you apply, warns you, and reduces the toothbrush vibrations automatically to help protect your gums. The toothbrush will make a

pulsing sound as a reminder to ease off the pressure. 7 out of 10 people found this feature helped them become better brushers.

SmartTimer and QuadPacer



2 minute SmarTimer and 30-second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

Designed around you



Slim, light-weight, and ergonomic design makes the toothbrush easy to hold and use so that you can get an amazing clean without any strain.

BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync

technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

Easy-Start



Our Easy-Start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

Long lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

Specifications

Power

Voltage: DC5V

Technical specifications

Operating time (full to empty): 14 days*****

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: Black

Service

Warranty: 2-year limited warranty

Ease of use

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

BrushSync Replacement Reminder:

Replacement reminder. Always know when to, replace brush heads, Reminder icon lights up, to always ensure best results

Items included

Brush heads: 1 W2 Optimal Black

Charger: 1

Handle: 1 3100 Series

Cleaning performance

Performance: Removes up to 3x more plaque*

Timer: Quadpacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday clean

Smart sensor technology

BrushSync Replacement Reminder: Always know when to, replace brush heads

Pressure sensor: Alerts when brushing too hard

