

# PHILIPS

## sonicare

Sonic electric  
toothbrush - Dispense

FlexCare Platinum

3 modes, 3 intensities

2 brush heads

With pressure sensor

UV brush head sanitizer



HX9182/32



## Their techniques may vary, but their results won't Our first brush with Adaptive Cleaning Technology

At the core of our Adaptive Cleaning Technology are two components - the innovative AdaptiveClean brush head and intuitive pressure sensor. Both features help your patients brush smarter, not harder.

### Proven to improve oral health

- Quadpacer
- Smarttimer
- UV sanitizer

### Designed around you

- InterCare brush head

### Guides you to follow dental recommendations

- Intuitive pressure sensor

### Provides a superior clean

- AdaptiveClean brush head
- Customizable brushing experiences
- Patented Sonicare toothbrush technology

### Offers a customized brushing experience

- New Deep Clean: Effortlessly tackle trouble spots


# Highlights

## AdaptiveClean brush head



Soft, flexible sides made of rubber allow the bristles to conform to the unique shape of patients' teeth and gums to provide up to 4x more surface contact\*. The brush head's flexibility helps the bristles gently track along the gumline and absorb excessive brushing pressure. Enhanced sweeping motion delivers a unique in-mouth feel and superior cleaning along the gumline and in hard-to-reach areas\*. Proven to remove 10x more plaque than a manual toothbrush in hard-to-reach areas\*.

## Intuitive pressure sensor




If too much pressure is applied, the handle gently pulsates to tell patients they are brushing too hard to guide proper technique. In a study of U.S. patients, 7 out of 10 aggressive brushers tested reduced their brushing pressure when using the pressure sensor.

## Deep Clean mode



Our all-new Deep Clean mode delivers an exceptional clean in hard-to-reach places. The unique motion of the brush head in this mode plus an extended three-minute brushing time makes it the ideal choice when you want to pay extra attention to trouble spots.

## Custom brushing experiences



Three intensity settings and three cleaning modes provide optimal comfort and performance for a more personalized brushing experience, encouraging better compliance. The three-minute Deep Clean mode delivers an exceptional clean in hard-to-reach places. It's the perfect mode for patients who need to pay extra attention to trouble spots.

## InterCare brush head



Extra-long, high-density bristles reach deeper between teeth to help remove significantly more plaque\* interproximally. Ideal for patients with bleeding gums or increased risk of gum inflammation.

## UV sanitizer




Kills up to 99% of bacteria and viruses\*. Can be available as an accessory in one of the dispensing SKUs.

## Smarttimer



Two-minute timer helps ensure dental professional recommended brushing time.

## Quadpacer



Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

# Specifications

## Design and finishing

**Color:** Platinum gray

## Power

**Voltage:** 110-220 V

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 3 weeks\*\*

**Battery type:** Lithium ION

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Battery indicator:** Light shows battery status

**Display:** Illuminated display

**Brush head system:** Easy click-on brush heads

**Handle:** Slim ergonomic design

**Brushing time:** Up to 3 weeks\*\*\*

## Items included

**Charger:** 1

**Handles:** 1 FlexCare Platinum

**Travel case:** 1

**Brush heads:** 1 AdaptiveClean, 1 InterCare standard

**UV sanitizer:** 1

## Cleaning performance

**Performance:** Removes up to 10x more plaque\*

**Health benefits:** Improves gum health in only two weeks

**Speed:** Up to 62000 brush movement/min

**Pressure sensor:** Vibrates handle to alert user

**Timer:** Quadpacer and SmarTimer

**Whitening benefits:** Helps naturally remove stains

## Modes

**3 intensities:** Low, Medium, High

**Clean:** For exceptional everyday clean

**White:** Removes surface stains

**Deep Clean:** For an invigorating deep clean



\* Than a DiamondClean brush head

\* Than a manual toothbrush

\* Data on file

\* After a single use of each toothbrush vs. Oral-B 3000/6000 CrossAction

\* HSV1, H.A.