

# PHILIPS

## sonicare

2 Series plaque control

1 mode

1 brush head



HX6211/18



## Removes up to 7X more plaque\* \*than a manual toothbrush

Offers you superior plaque removal at the push of a button

### **Guides you to follow dentist recommendations**

- SmarTimer helps you meet recommended brushing time

### **Offers a customized brushing experience**

- EasyStart program helps you adjust from a manual toothbrush

### **Provides a superior clean you can see and feel**

- Removes up to 7x more plaque than a manual toothbrush
- Angled brush head for better access to back teeth
- Better cleaning in hard-to-reach areas than a manual brush
- Sonicare dynamic cleaning action drives fluid between teeth

### **Recommended by dental professionals world wide**

- Better plaque removal to help reduce cavities
- Safe on orthodontics, dental work and veneers

# Highlights

## Helps reduce cavities



Provides better plaque removal and access to back teeth; brushing twice a day every day with this brush can help reduce cavities

## Up to 7x better plaque removal



Removes up to 7x more plaque than a manual toothbrush for a superior clean

## Sonic technology



Our sonic technology creates a unique dynamic cleaning action that drives fluid deep between teeth and along the gumline for a thorough-yet-gentle clean

## Safe and gentle




The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worry.

## Better access to back teeth




This plaque-removing toothbrush has a slim, angled brush head which provides better access to back teeth and hard-to-reach places.

## Better in hard-to-reach areas



The slim, angled neck and easy-grip handle of this toothbrush helps with better cleaning in hard-to-reach areas than a manual toothbrush

## EasyStart



The EasyStart feature adjusts the power level of the toothbrush over the first 14 uses to help you get used to the Sonicare experience, streamlining the transition from a manual toothbrush

## SmarTimer



This Philips Sonicare power toothbrush features a special SmarTimer to help you meet the dentist-recommended brushing time of two minutes

# Specifications

## Power

**Voltage:** 110–220 V, 110–220 V

## Technical specifications

**Operating time (full to empty):** Up to 2 weeks, Up to 2 weeks

**Battery:** Rechargeable, Rechargeable

**Battery type:** NiMH, NiMH

## Design and finishing

**Color:** Cambridge Blue Green, Cambridge Blue Green

## Service

**Warranty:** 2-year limited warranty, 2-year limited warranty

## Ease of use

**Battery indicator:** Light shows battery status, Light shows battery status

**Brush head system:** Easy click-on brush heads, Easy click-on brush heads

**Handle:** Slim ergonomic design, Slim ergonomic design

**Brushing time:** Up to 2 weeks, Up to 2 weeks

## Items included

**Handles:** 1 2 Series plaque control

**Charger:** 1, 1

**Brush heads:** 1 ProResults plaque control, 1 ProResults plaque control

## Cleaning performance

**Performance:** Removes up to 7x more plaque\*, Removes up to 7x more plaque\*

**Health benefits:** Helps reduce cavities, Helps reduce cavities

**Speed:** Up to 62000 brush movement/min, Up to 62000 brush movement/min

**Timer:** SmarTimer, SmarTimer

## Modes

**Clean:** For exceptional everyday clean, For exceptional everyday clean

