

# PHILIPS

## AVENT

DECT Baby Monitor

Audio Monitors

100% private connection

Starry night light projector

Talkback and MP3 lullabies

Vibration alert

 Product Picture

SCD580/01



## The most reliable connection to your baby

Reassurance and well-being for you and your baby

Our Philips Avent SCD580/01 Baby Monitor delivers complete comfort and reassurance for you and your baby. Providing the most reliable connection combined with a variety of soothing features. Your baby will love the night light projector!

### Freedom and flexibility

- Range up to 330 meters\*
- Superior operating time up to 18 hours
- Vibration alert
- Instant dim mode for a peaceful night

### Comforting and soothing your baby

- Stay connected with talk back
- Night light + MP3 plug and play
- Calm your baby with the starry night light projector

### Reliable, secure connection & constant reassurance

- Always know when your monitor is in range and connected
- Perfect clear sound due to DECT Technology
- Energy saving Smart ECO mode
- Sound level lights indicate if the baby makes noise
- Monitor the humidity and temperature in your baby's room

# Highlights

## Starry night light projector

Your baby may still be restless when you first put them to bed. Help soothe and calm your baby to sleep with the starry night light projection which can be activated from the parent or baby unit. (Remote activation not available in US and Canada).

## Night light with MP3 plug

Gently calm your baby with your own selected song and a warm tranquil glow of the night light. With MP3 plug and play you can play your own selected songs. Your personal music in the baby room helping your baby drift to sleep in no time. (Remote activation not available in US and Canada).

## Zero interference with DECT



DECT technology guarantees zero interference from any other transmitting product like other baby monitors, cordless phones and cell phones. Data encryption provides a secure and private connection, so you can be sure that you are the only one that can hear your baby.

## Perfect clear sound



Hear every tiny giggle, gurgle, and hic-up with perfect clarity. DECT (Digital Enhanced Cordless Telecommunications) technology provides high quality, crystal clear sound so you can hear your baby anytime.

## Smart ECO mode




The unique Smart ECO mode automatically minimizes your transmission power and increases your battery lifetime. The closer you are to your baby, the less power is needed for a perfect connection (not available in US and Canada).

## Sound activation 5 LED lights



Lights are activated to indicate the level of sound in your baby's room. Even when the Parent Unit is muted.

## Always connected



The parent unit will alert you when your monitor is out of range or when the power is low, helping you make sure that you are always connected to your baby.

## Temperature & humidity sensor

Combined humidity and temperature sensor allows you to monitor the climate in your baby's room. Your baby can't regulate its body

temperature as well as you can and a slight change in temperature, and too much or too little humidity can make your baby restless. A specially combined humidity and temperature sensor will immediately alert you, through the digital display on the parent unit, if the climate in your baby's room changes. Monitor the complete climate in baby's room.

## 330m Range



Indoor range up to 50 meters\*. Outdoor range up to 330 meters\*.

## Operating time up to 18 hours



The small, rechargeable parent unit will give you the freedom of cordless streaming for up to 18 hours before re-charging is required.

## Vibration alert

Personalize your alert settings and get notified by a silent and subtle vibration of the parent unit.

## Peaceful night mode

Press the night dim button to instantly activate the peaceful night mode. In this mode the light intensity of the lights, display as well as the volume of the alerts are reduced. Easily press the night dim button again in order to switch back to normal operating mode.

# Specifications

## Development stages

Stage: Pregnancy, 0 - 6 months

## Features

Comforting night light  
Indoor range up to: 50m  
Linked and in range indication  
Outdoor range up to: 330 meters  
Talkback function  
LCD Screen  
Soothing lullabies

## Power

Operating time on battery hour: 18 hours  
Charging time: 10 hour(s)  
Power Supply: 120 V (US), 220 - 240 V

## Tuner/Reception/Transmission

Automatic channel selection  
Frequency band: DECT

## Convenience

Auto out-of-range warning  
Battery low indication  
Battery charging indication  
Volume control  
Sound-level lights

## Convenience

Sensitivity control

## Technical specifications

Storage temperature range: 10 - 40 °C  
Operating temperature range °C: 10 - 40 °C

## Logistic data

F-box dimensions (W x H x D): 220 x 156 x 94 mm

## Features

Charging station  
DECT Technology  
Energy saving Smart ECO mode: ECO mode, Yes

Peaceful night mode and vibration alert  
Sound activation lights: 5 LED lights  
Starry night light projector  
Temperature/humidity sensor

## Consumer Trade Item

Height: 15.6  
Width: 22  
Length: 9.6  
Net Weight: .877  
Gross Weight: 1.045  
GTIN: 08710103610908  
Country of origin: CN  
Harmonized system code: 851762

## Outer Carton

Height: 18.3  
Width: 23.5  
Length: 30.1  
Gross Weight: 3.445  
GTIN: 18710103610943

