

PHILIPS

sonicare

Sonic electric
toothbrush - Dispense

FlexCare Platinum

3 modes, 3 intensities

1 brush head

With pressure sensor

UV brush head sanitizer

HX9181/25



Their techniques may vary, but their results won't Our first brush with Adaptive Cleaning Technology

At the core of this technology are two components: the AdaptiveClean brush head and an intuitive pressure sensor. Both features help your patients to brush smarter.

Proven to improve oral health

- Quadpacer
- Smarttimer
- UV sanitizer

Guides you to follow dental recommendations

- Intuitive pressure sensor

Provides a superior clean

- AdaptiveClean brush head
- Patented Sonicare toothbrush technology
- Customizable brushing experiences

Highlights

Sonic technology




Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

AdaptiveClean brush head




Soft, flexible sides made of rubber allow the bristles to conform to the unique shape of patients' teeth and gums to provide up to 4x more surface contact*. The brush head's flexibility helps the bristles gently track along the gumline and absorb excessive brushing pressure. Enhanced sweeping motion delivers a unique in-mouth feel and superior cleaning along the gumline and in hard-to-reach areas*. Proven to remove 10x more plaque than a manual toothbrush in hard-to-reach areas*.

Intuitive pressure sensor



If too much pressure is applied, the handle gently pulsates to tell patients they are brushing too hard to guide proper technique. In a study of U.S. patients, 7 out of 10 aggressive brushers tested reduced their brushing pressure when using the pressure sensor.

Custom brushing experiences




Three intensity settings and three cleaning modes provide optimal comfort and performance for a more personalized brushing experience, encouraging better compliance. The three-minute Deep Clean mode delivers an exceptional clean in hard-to-reach places. It's the perfect mode for patients who need to pay extra attention to trouble spots.

UV sanitizer



Kills up to 99% of bacteria and viruses*. Can be available as an accessory in one of the dispensing SKUs.

Quadpacer



Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

Smarttimer



Two-minute timer helps ensure dental professional recommended brushing time.

Specifications

Design and finishing

Color: Platinum gray

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 3 weeks**

Battery: Rechargeable

Battery type: Lithium ION

Service

Warranty: 2-year limited warranty

Ease of use

Display: Illuminated display

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design

Brushing time: Up to 3 weeks***

Items included

Handles: 1 FlexCare Platinum

Charger: 1

Travel case: 1

Brush heads: 1 InterCare standard

UV sanitizer: 1

Cleaning performance

Performance: Removes up to 7x more plaque*

Health benefits: Improves gum health in only two weeks

Speed: Up to 62000 brush movement/min

Pressure sensor: Vibrates handle to alert user

Timer: Quadpacer and SmartTimer

Whitening benefits: Helps naturally remove stains

Modes

Clean: For exceptional everyday clean

3 intensities: Low, Medium, High

White: Removes surface stains

Deep Clean: For an invigorating deep clean



* Than a DiamondClean brush head

* Than a manual toothbrush

* Data on file

** based on two periods of two-minute brushings per day