

# PHILIPS

## sonicare

Sonic electric  
toothbrush - Dispense

FlexCare Platinum

3 modes, 3 intensities

1 brush head

With pressure sensor

UV brush head sanitizer

HX9181/25



## Their techniques may vary, but their results won't

### Our first brush with Adaptive Cleaning Technology

At the core of this technology are two components: the AdaptiveClean brush head and an intuitive pressure sensor. Both features help your patients to brush smarter.

#### **Proven to improve oral health**

- Quadpacer
- Smarttimer
- UV sanitizer

#### **Guides you to follow dental recommendations**

- Intuitive pressure sensor

#### **Provides a superior clean**

- AdaptiveClean brush head
- Patented Sonicare toothbrush technology
- Customizable brushing experiences

# Highlights

## Sonic technology




Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

## AdaptiveClean brush head




Soft, flexible sides made of rubber allow the bristles to conform to the unique shape of patients' teeth and gums to provide up to 4x more surface contact\*. The brush head's flexibility helps the bristles gently track along the gumline and absorb excessive brushing pressure. Enhanced sweeping motion delivers a unique in-mouth feel and superior cleaning along the gumline and in hard-to-reach areas\*. Proven to remove 10x more plaque than a manual toothbrush in hard-to-reach areas\*.

## Intuitive pressure sensor



If too much pressure is applied, the handle gently pulsates to tell patients they are brushing too hard to guide proper technique. In a study of U.S. patients, 7 out of 10 aggressive brushers tested reduced their brushing pressure when using the pressure sensor.

## Custom brushing experiences




Three intensity settings and three cleaning modes provide optimal comfort and performance for a more personalized brushing experience, encouraging better compliance. The three-minute Deep Clean mode delivers an exceptional clean in hard-to-reach places. It's the perfect mode for patients who need to pay extra attention to trouble spots.

## UV sanitizer



Kills up to 99% of bacteria and viruses\*. Can be available as an accessory in one of the dispensing SKUs.

## Quadpacer



Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

## Smarttimer



Two-minute timer helps ensure dental professional recommended brushing time.

# Specifications

## Power

**Voltage:** 110–220 V

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 3 weeks\*\*

**Battery type:** Lithium ION

## Design and finishing

**Color:** Platinum gray

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Display:** Illuminated display

**Battery indicator:** Light shows battery status

**Brush head system:** Easy click-on brush

heads

**Handle:** Slim ergonomic design

**Brushing time:** Up to 3 weeks\*\*\*

## Items included

**Charger:** 1

**Handles:** 1 FlexCare Platinum

**Travel case:** 1

**Brush heads:** 1 InterCare standard

**UV sanitizer:** 1

## Cleaning performance

**Performance:** Removes up to 7x more plaque\*

**Health benefits:** Improves gum health in only two weeks

**Speed:** Up to 62000 brush movement/min

**Pressure sensor:** Vibrates handle to alert user

**Timer:** Quadpacer and SmarTimer

**Whitening benefits:** Helps naturally remove stains

## Modes

**Clean:** For exceptional everyday clean

**3 intensities:** Low, Medium, High

**White:** Removes surface stains

**Deep Clean:** For an invigorating deep clean

