

# PHILIPS

## sonicare

Sonic electric  
toothbrush

HealthyWhite

2 modes

1 brush head

HX6711/07



## Whiter, healthier teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. HealthyWhite is proven to remove everyday stains in just two weeks, with the regular use of the Clean & White Mode.

### Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology

### Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

### Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

### Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2 minute mode

# Highlights

## Sonic technology




Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## SmarterTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.\*

## Clean Mode

For exceptional everyday clean. Helps maintain whitening treatments.



# Specifications

## Design and finishing

**Color:** Glacier Green

## Power

**Voltage:** Multi-voltage charger

## Technical specifications

**Battery:** Rechargeable

## Features

**ProResults brush head**

**Smartimer:** 2-minute timer helps ensure recommended brushing

**Quadpacer:** Encourages attentive brushing of each quadrant of the mouth

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Handle:** Slim ergonomic design and rubber grip for easy handling.

**Brush head system:** Easy snap-on brush heads for optimal hygiene

**Battery indicator:** Illuminated icon blinks when its time to recharge

**Brushing time:** 42 2-minute brushings OR 3 weeks

## Items included

**Charger:** 1 multi-voltage charger

**Charger cover:** 1 charger cover with brush head holder

## Cleaning performance

**Brush heads:** 1 ProResults standard brush head for an All-Around Clean

**Performance:** Removes more plaque than Sonicare Essence/Sonicare CleanCare

**Health benefits:** Helps improve gum health

**Speed:** Up to 62000 brush movement/min

**Whitening benefits:** Can whiten teeth up to two shades in two weeks

**Timer:** Quadpacer and Smartimer

**Modes:** 2 Modes: Clean and Clean and White

