



Philips Sonicare
AirFloss Pro/Ultra -
Dispense

2 nozzles

2 oz. BreathRx mouth rinse



HX8382/14

A triple burst of better gingival health

Efficacy meets compliance for infrequent flossers

Philips Sonicare AirFloss Pro/Ultra**** is clinically proven as effective as floss in improving gingival health.* In fact, up to 97% of inconsistent flossers who tried AirFloss Pro/Ultra showed improved gum health.*

Superior results for patients' oral care needs

- Exclusive bonus items only available to dental professionals

Improves patients' oral health

- Clinically proven as effective as floss for gum health**
- Helps prevent cavities between teeth
- Up to 99.9% plaque removal***

Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

The easiest way to effectively clean between teeth

- Takes just 60 seconds: Point, Press, Clean!
- Fresh and clean with mouthwash
- An easy way to start a healthy habit

PHILIPS

sonicare


Highlights

Up to 99.9% plaque removal



AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.**

Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Cleans in just 60 seconds



With the AirFloss Pro/Ultra, patients can clean their entire mouth in less than 60 seconds a day. Patients simply select their desired burst frequency (single, double or triple) and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Fresh and clean with mouthwash



Patients can use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Helps prevent cavities between



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps


prevent cavities from forming in the spaces between patients' teeth.

High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your patient's liking with single, double, or triple bursts for each press of the activation button.

Bonus items not sold at retail

Now you can offer these exclusive Philips Sonicare AirFloss Pro/Ultra bonus items your patients won't find anywhere else. Designed around patient needs for dental offices like yours, the package includes: 1 extra nozzle and travel bag.

Specifications

Design and finishing

- Color: Pink

Power

- Voltage: Multi-voltage charger

Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

Service

- Warranty: 2-year limited warranty

Ease of use

- Battery indicator: Handle LEDs show charge level
- Handle: Slim ergonomic design
- Battery Life: 1-3 weeks*****
- Nozzle attachment: Easily snaps on and off

Items included

- Charger: 1
- Antibacterial mouth rinse: 2oz
- Nozzle holder: 1
- Travel bag: 1
- AirFloss Pro/Ultra handle: 1
- AirFloss Pro/Ultra nozzle: 2

Cleaning performance

- Performance: Removes up to 99.9% plaque**
- Health benefits: Improve gum health in 2 weeks*, Cleans entire mouth in 60 sec
- Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction
- Customizable settings: Auto-burst or manual burst, Single/double/triple burst***
- For best results: Change nozzle every 6 months

Ingredients

- Antibacterial mouth rinse: 0.075% CPC, Zytex



Issue date 2022-06-05

Version: 2.2.1

EAN: 00 07502 00584 85

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* When used in conjunction with a manual toothbrush and antimicrobial rinse in patients with mild to moderate gingivitis. AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. For more information, please visit www.philips.com/airfloss/faq.

*** AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.
***** depending on burst setting used