

PHILIPS

sonicare

2 Series plaque control

2 Series plaque control

1 mode

2 brush heads

2 handles



HX6252/71



Removes up to 6X more plaque* *than a manual toothbrush

Offers you superior plaque removal at the push of a button

Guides you to follow dentist recommendations

- Quadpacer interval timer encourages thorough brushing
- SmarTimer helps you meet recommended brushing time

Offers a customized brushing experience

- EasyStart program helps you adjust from a manual toothbrush

Provides a superior clean you can see and feel

- Removes up to 7x more plaque than a manual toothbrush
- Angled brush head for better access to back teeth
- Better cleaning in hard-to-reach areas than a manual brush
- Sonicare dynamic cleaning action drives fluid between teeth

Recommended by dental professionals world wide

- Better plaque removal to help reduce cavities
- Safe on orthodontics, dental work and veneers

Highlights

Helps reduce cavities



Provides better plaque removal and access to back teeth; brushing twice a day every day with this brush can help reduce cavities

Up to 7x better plaque removal



Removes up to 7x more plaque than a manual toothbrush for a superior clean

Sonic technology



Our sonic technology creates a unique dynamic cleaning action that drives fluid deep between teeth and along the gumline for a thorough-yet-gentle clean

Safe and gentle




The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worry.

Better access to back teeth




This plaque-removing toothbrush has a slim, angled brush head which provides better access to back teeth and hard-to-reach places.

Better in hard-to-reach areas



The slim, angled neck and easy-grip handle of this toothbrush helps with better cleaning in hard-to-reach areas than a manual toothbrush

EasyStart



The EasyStart feature adjusts the power level of the toothbrush over the first 14 uses to help you get used to the Sonicare experience, streamlining the transition from a manual toothbrush

SmarTimer



This Philips Sonicare power toothbrush features a special SmarTimer to help you meet the dentist-recommended brushing time of two minutes

Quadpacer



Features a Quadpacer interval timer to help encourage thorough brushing of the entire mouth. The toothbrush beeps every 30 seconds, signifying it's time to move to another area of the mouth.

Specifications

Power

Voltage: 110–220 V, 110–220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks, Up to 2 weeks

Battery: Rechargeable, Rechargeable

Battery type: NiMH, NiMH

Design and finishing

Color: Sky Blue and Ultra Coral, Sky Blue and Ultra Coral

Service

Warranty: 2-year limited warranty, 2-year limited warranty

Ease of use

Battery indicator: Light shows battery status, Light shows battery status

Brush head system: Easy click-on brush heads, Easy click-on brush heads

Handle: Slim ergonomic design, Slim ergonomic design

Brushing time: Up to 2 weeks, Up to 2 weeks

Items included

Charger: 2, 2

Handles: 2 2 Series plaque control

Travel case: 2, 2

Brush heads: 2 ProResults plaque control, 2 ProResults plaque control

Cleaning performance

Performance: Removes up to 6X more plaque*, Removes up to 7x more plaque*

Health benefits: Helps reduce cavities, Helps reduce cavities

Speed: Up to 31000 brush movement/min, Up to 31000 brush movement/min

Timer: Quadpacer and SmarTimer, Quadpacer and SmarTimer

Modes

Clean: For exceptional everyday clean, For exceptional everyday clean

