

# PHILIPS

## sonicare

Sonic electric  
toothbrush

HealthyWhite

2 modes

1 brush head

HX6711/09



## Whiter, healthier teeth

### The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Sonicare HealthyWhite electric toothbrush HX6711/02 is proven to remove everyday stains in just two weeks, with the use of the Clean & White Mode.

#### **Proven to improve oral health**

- Philips Sonicare toothbrush helps whiten teeth

#### **Guides you to follow dental recommendations**

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

#### **Provides a superior clean**

- Sonicare dynamic cleaning action drives fluid between teeth
- Removes up to 2x more plaque than a manual toothbrush

#### **Offers a customized brushing experience**

- Easy-start program to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains

# Highlights

## Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

## Helps whiten teeth



This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

## Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.\*

## SmarterTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.



# Specifications

## Design and finishing

**Color:** Glacier Green

## Power

**Voltage:** 110-220 V

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 3 weeks\*\*

**Battery type:** Lithium ION

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Handle:** Slim ergonomic design, Rubber grip for easy handling

**Brush head system:** Easy click-on brush heads

**Battery indicator:** Light shows battery status

**Brushing time:** Up to 3 weeks\*\*

## Items included

**Brush heads:** 1 ProResults standard

**Charger:** 1

**Handles:** 1 HealthyWhite

## Cleaning performance

**Performance:** Removes up to 90% of stains\*

**Speed:** Up to 62000 brush movement/min

**Health benefits:** Helps improve gum health

**Whitening benefits:** Whitens teeth up to 2 shades

**Timer:** Quadpacer and SmarTimer

## Modes

**Clean:** For exceptional everyday clean

**Clean & White:** 2 min Clean + 30 sec White

