

# PHILIPS

## sonicare

2 Series plaque  
control

1 mode

2 brush heads

2 handles



HX6212/19



## Removes up to 2X more plaque\*

\*than a manual toothbrush

Offers you superior plaque removal at the push of a button

### Proven to improve oral health

- Safe and gentle for your oral health needs

### Designed around you

- Easy start program builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing

### Provides a superior clean

- Removes up to 2x more plaque than a manual toothbrush
- Better plaque removal to help reduce cavities

### Optimize your brushing

- Philips Sonicare's advanced sonic technology
- Angled brush head for better access to back teeth

# Highlights

## Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

## Helps reduce cavities




Provides better plaque removal and access to back teeth; brushing twice a day every day with this brush can help reduce cavities

## The only true sonic technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

## Better access to back teeth




This plaque-removing toothbrush has a slim, angled brush head which provides better access to back teeth and hard-to-reach places.

## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

## Safe and gentle



Whatever your specific needs, you can be sure of a safe brushing experience with the 2 Series plaque control. Our sonic technology is suitable for use with braces, fillings, crowns and veneers, and can even be used in the daily treatment of gum disease.

# Specifications

## Power

**Voltage:** 110–220 V

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 2 weeks

**Battery type:** NiMH

## Design and finishing

**Color:** Cambridge Blue Green, Black

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Battery indicator:** Light shows battery status

**Brush head system:** Easy click-on brush heads

**Handle:** Slim ergonomic design

**Brushing time:** Up to 2 weeks

## Items included

**Charger:** 1

**Brush heads:** 2 C1 SimplyClean standard

**Handle:** 2 2-series plaque control

## Cleaning performance

**Performance:** Up to 2X more plaque removal\*

**Health benefits:** Helps reduce cavities

**Speed:** Up to 31000 brush movement/min

**Timer:** Quadpacer and SmarTimer

## Modes

**Clean:** For exceptional everyday clean

