

IHE Integration Statement forNotify



| Vendor | Product Name | Version | Date |
|---------|--------------|---------|------------|
| Forcare | forNotify | 2022-2 | 2022-05-20 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

| Integration Profile | Actor | Option |
|--|--|------------------------|
| Consistent Time (CT) | Time Client | |
| Audit Trail and Node Authentication (ATNA) | Secure Application (ITI-19, ITI-1, ITI-20) | |
| Cross-Enterprise User Assertion (XUA) | X-Service-Provider (ITI-40) | Subject-Role |
| | | PurposeOfUse |
| | X-Service-User (ITI-40) | Subject-Role |
| | | PurposeOfUse |
| Internet User Authorization (IUA) | Resource Server (ITI-72) | SAML Token |
| Document Metadata Subscription (DSUB) | Document Metadata Notification Recipient (ITI-53) | |
| | Document Metadata Notification Subscriber (ITI-52) | |
| Cross-Enterprise Document Sharing (XDS.b) | Document Consumer (ITI-18, ITI-43) | |
| Cross-Enterprise Document Workflow (XDW) | Content Consumer | Document Import Option |

Internet address for vendor's IHE information www.forcare.com

Links to Standards Conformance Statements for the Implementation

| | |
|-------|--|
| HL7 | www.forcare.com |
| DICOM | www.forcare.com |

Links to general information on IHE

In North America: www.ihe.net In Europe: www.ihe-europe.net In Japan: www.ihe-j.org

© 2006-2022 Forcare Holding B.V. All rights reserved.