

# PHILIPS

## sonicare

Sonic electric  
toothbrush

FlexCare Platinum

3 modes, 3 intensities

2 brush heads

With pressure sensor

HX9110/02



## Exceptional clean between teeth

Removes up to 7x more plaque between teeth\*

With nine unique brushing options, an intuitive pressure sensor and advanced brush head technology, the Philips Sonicare FlexCare Platinum toothbrush goes beyond traditional cleaning to deliver exceptional plaque removal and improved gum health.

### Proven to improve oral health

- Improves gum health in just 2 weeks
- Safe and gentle for your oral health needs

### Designed around you

- Timers encourage thorough 2 minutes of brushing
- Long battery life and travel case
- Alerts you if you're pressing too hard

### Provides a superior clean

- Our unique technology brings you a powerful yet gentle clean
- Remove up to 7x more plaque with our InterCare brush head\*

### Optimize your brushing

- Choose from 3 modes and 3 intensity settings


# Highlights

## Remove up to 7x more plaque




Click on our InterCare brush head to improve gum health in just 2 weeks. Extra-long bristles help remove more plaque in hard-to-reach spots and between teeth for healthy gums.

## Look after those gums




With optimal cleaning from your FlexCare, gums will become healthier in 2 weeks. Removing up to 7x more plaque along the gum line than a manual toothbrush, you'll get your healthiest smile.

## 3 modes and 3 intensity settings



With FlexCare Platinum, you experience a truly deep clean. Our 3 intensity settings let you boost your clean, while 3 modes meet your specific needs: Clean mode – for superior daily cleaning, White – to remove surface stains, and Clean for outstanding daily cleaning.

## Too much pressure?



You might not notice if you're brushing too hard, but your FlexCare Platinum will. If you use too much pressure, your handle will gently pulsate. This is a reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

## Philips Sonicare technology



Powerful sonic vibrations whip up your toothpaste into plaque fighting bubbles, and drive them deep between your teeth and along your gum line. At the same time your teeth will experience 62,000 gentle but effective brush movements. You'll get a month's worth of manual brushing in just 2 minutes.

## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of

time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Brush with confidence



Whatever your specific needs, you can be sure of a safe brushing experience with FlexCare Platinum. Our sonic technology is suitable for use with braces, fillings, crowns and veneers, and can even be used in the daily treatment of gum disease.

## Designed to fit your lifestyle



Our premium travel case lets you store your toothbrush hygienically, while our compact travel charger keeps you topped up when you're on the go. You'll enjoy two weeks of regular use from a single full charge, so being on the move means you needn't miss out on that fresh feeling.

# Specifications

## Power

**Voltage:** 110–220 V

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 3 weeks\*\*

**Battery type:** Lithium ION

## Design and finishing

**Color:** Platinum gray

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Display:** Illuminated display

**Battery indicator:** Light shows battery status

**Brush head system:** Easy click-on brush heads

**Handle:** Slim ergonomic design

**Brushing time:** Up to 3 weeks\*\*

## Items included

**Charger:** 1

**Travel case:** 1

**Brush heads:** 2 InterCare standard

**Handle:** 1 FlexCare Platinum

## Cleaning performance

**Performance:** Removes up to 7x more plaque\*

**Health benefits:** Improves gum health in only two weeks

**Speed:** Up to 62000 brush movement/min

**Pressure sensor:** Vibrates handle to alert user

**Timer:** Quadpacer and SmarTimer

**Whitening benefits:** Helps naturally remove stains

## Modes

**Clean:** For exceptional everyday clean

**3 intensities:** Low, Medium, High

**White:** Removes surface stains

**Gum Care:** Gently massages gums

