



Philips Sonicare FlexCare  
Sonic electric toothbrush

### 3 modes

2 brush heads  
UV brush head sanitizer



HX6982/10

## Advanced Cleaning

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare that adapts to your oral care needs.

#### Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology
- Improves gum health in only two weeks

#### Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

#### Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- UV brush head Sanitizer - Cleaner brush, better oral health

#### Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience
- Three different modes for a better cleaning experience

**PHILIPS**  
**sonicare**

# Highlights

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck




The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## UV brush head Sanitizer

All-in-one sanitizing, charging and storage. The UV brush head Sanitizer helps keep your brush head free of bacteria. It kills up to 99% of bacteria and viruses.\*

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## Improves gum health




This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

## SmarTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer




30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Three different cleaning modes



The three modes include a clean mode for an outstanding clean in 2 minutes. A sensitive mode for gentle yet thorough cleaning of the gums and a massage mode to stimulate the gum.



Sonic electric toothbrush  
3 modes 2 brush heads, UV brush head sanitizer

HX6982/10

# Specifications

## Design and finishing

- Color: White & green

## Power

- Voltage: 110-220 V

## Technical specifications

- Operating time (full to empty): Up to 2 weeks
- Battery: Rechargeable
- Battery type: Lithium ION

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Brush head system: Easy click-on brush heads
- Handle: Slim ergonomic design, Rubber grip for easy handling
- Battery indicator: Light shows battery status
- Brushing time: Up to 2 weeks

## Items included

- Travel case: 1 soft
- Charger: 1
- Handles: 1 FlexCare
- Brush heads: 1 ProResults standard, 1 ProResults compact
- UV sanitizer

## Cleaning performance

- Speed: Up to 62000 brush movement/min
- Performance: Removes up to 83% of plaque in hard-to-reach areas
- Health benefits: Helps improve gum health
- Whitening benefits: Whitens teeth up to 2 shades
- Timer: Quadpacer and SmarTimer

## Modes

- Clean: For exceptional everyday clean
- Massage: Invigorating massage
- Sensitive: Gentle teeth and gum cleaning
- 2 Routines: Go Care and Max Care



Issue date 2023-01-14

Version: 3.3.2

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* E coli, S. mutans and Herpes simplex