

PHILIPS

AVENT

Single electric
breast pump

Ultra Comfort

Includes 4oz Natural bottle



SCF332/61



More comfort, more milk*

Clinically proven comfort*

Sit comfortably without leaning forward as our soft massage cushion gently stimulates your milk flow. Use our quiet pump anytime, anywhere, even with batteries. It's easy to set up, personalize, use and clean.

Effective expression

- 3 customizable settings at your fingertips

Easy to set up, express, clean and travel with

- Easy to use anywhere
- Intuitive design, easy to assemble and clean
- No spills or messy transfers

Express in comfort

- Clinically proven less leaning forward*
- Designed to let you express quietly
- Soft cushion with massaging petals

The support in your breastfeeding journey

- Compatible with our natural bottles, for easy feeding
- Easy access to warranty and support

Highlights

Relaxing* expression position



Research* suggests that when you're relaxed, your milk flows more easily. Thanks to its unique design, you won't need to learn forward when you pump. With this breast pump, you can comfortably sit up while your milk flows directly from your breast into the bottle.

Soft massage cushion



The massage cushion is made of a soft, velvety material to give you a feeling of warmth and helps stimulate your milk flow. With nature in mind, the iconic petal cushion is designed to mimic baby suckling, which is intended to gently help stimulate let down.

Quiet pump

Designed for discreet expression, anytime.

Control breast milk flow

With the press of a button, you can start pumping in the gentle stimulation mode to activate your let-down. To begin breast milk collection, choose one of the three flow

settings for your preferred suction speed. To increase the suction, press a button again to adjust.

Compact and portable

Our expression kit, which includes the funnel, cushion and bottle, is easy to hold in one hand and position on your breast. Place the lightweight base unit within easy reach for full control when pumping. It's also very convenient when you go out. Just wrap the tube around the base unit. If you plan to travel abroad, read the specifications tab to check the voltage before plugging it in or use the batteries for easy back-up.

Intuitive assembly



The breast pump features an intuitive design that is easy to clean. And since it only has a few parts, you can quickly assemble it and get started. Before each use, we recommend you wash the parts in warm, soapy water and rinse. Then sanitize with a Philips Avent steam sterilizer, or boil them for five minutes for best results. Use a damp cloth to clean the motor unit, cap, silicone tube and adapter.

No milk spillage in the tubes



Thanks to its closed system, the breast pump collects your milk directly into the storage container. This way your milk is collected while the suction tubes remain clean.

Handy accessories



To help you and baby on your breastfeeding journey, you'll find a 4oz Philips Natural bottle included with your pump. The wide, breast-shaped teat promotes natural latch-on similar to the breast to help make it easier for your baby to combine breast with bottle feeding. The pump also fits our breast feeding accessories, such as milk storage cups. See our Accessories tab for more options and replacement parts.

Warranty coverage

Your breast feeding journey is an ongoing experience. That's why we offer specialized product warranty coverage. Register your breast pump for product coverage as well as special offers and easy access to product support. Click the Support tab to see eligible details in your country.

Specifications

Design

Bottle design: Ergonomic shape, Wide neck
Breast pump design: Compact design

Country of origin

England

Material

Bottle: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

Breast pump: BPA free* (food contact parts only)

What is included

Sealing disc for milk storage: 1 pcs

Natural bottle 4oz: 1 pcs

Breast pump body: 1 pcs

Base unit incl. tubing: 1 pcs

Travel cover: 1 pcs

Standard size cushion(19.5 mm): 1 pcs (a larger size cushion is available and sold separately)

Disposable breast pads: 4 pcs

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Settings: 1 Stimulation mode, 3 Expression settings

Development stages

Stage: 0 – 6 months



* Clinically proven comfort: In a test carried out on 110 mothers in the USA, UK, China and Russia in March 2016, mothers awarded an average of 8.6/10 score for Philips Avent performance on comfort.

* More milk: Independent research has shown that there may be a link between higher comfort levels and increased milk production. See www.philips.com/AVENT

* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011