

Philips Philips Avent
Manual breast pump
diaphragm

White
Breast pumps

CP1970/01



Creates the vacuum needed for expression

Check for compatibility below

This diaphragm creates the vacuum needed for expressing breast milk. For use with manual breast pumps.

PHILIPS

Specifications

Replaceable parts

• Fit product types: SCF430/01 SCF430/03 SCF430/

10 SCF430/13 SCF430/15 SCF430/16 SCF430/20
SCF430/30 SCF441/01



Issue date 2024-04-13

Version: 5.5.1

© 2024 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com