

PHILIPS

AVENT



Electric breast pump Premium Plus

Double

Premium plus

Rechargeable battery

Timer display



SCF394/11

Inspired by baby. Effective for mum

Natural motion technology for quicker milk flow*

Enter a new era of expression with the Philips Avent Electric Breast Pump. Perfect balance of suction & nipple stimulation, inspired by baby's drinking rhythm and soft cushion gently adapts to your nipple size & shape. See more below.

Other benefits

- Choose the right nipple flow for your baby

Personalized and easy to use

- Minimal parts and intuitive setup
- Pause/Play function
- Memory function
- Quiet motor for a discreet experience anywhere
- Breast pump bag and pouch
- Personalized experience, 8 + 16 setting levels

Gentle and comfortable


- One size soft and adaptive silicone cushion
- Express without leaning forward

Quick expressing, more milk in less time*

Rechargeable battery and timer display


Highlights

Quick milk flow*



Express more milk in less time* with a cushion that stimulates the breast to express milk just like baby. It seamlessly adjusts from stimulation mode to expression mode¹, and applies just the right amount of nipple stimulation and suction for maximum milk flow. Based on milk flow initiation time (time to Milk Ejection Reflex - MER) results.*

Adapts to your nipple



One size fits all. Because we all come in different shapes and sizes, the silicone cushion gently flexes and adapts to fit your nipple. It fits 99.98% of nipple sizes* (up to 30mm).

Finding the right nipple is important



If your newborn baby is consistently not taking enough milk throughout the feeding sessions or has complications in getting milk, switch to a nipple with a higher flow rate. If persistent feeding issues occur, consult a healthcare professional.

Find the right nipple flow



Every baby feeds differently, and develops at their own pace. We've designed a range of flow rates so you can find the perfect one for your baby and personalize your bottle. All of the Natural Response Nipples are made of soft silicone.

Sit comfortably



Feel relaxed as you express thanks to a design that lets you sit upright rather than having to lean forward. No need to lean forward, clinically proven comfortable expressing position*. Based on product clinical trial results with 20 participants (2019); 90% of participants find the expressing position comfortable (single electric); 95% of participants find the expressing position comfortable double electric).

Personalize to your needs



Finely tune every session to your needs with a wide range of stimulation and expression settings. Our breast pump offers 8 stimulation and 16 expression levels for a personalized experience.

Express anywhere discreetly



Express whenever and wherever you like thanks to our quiet motor with a slimline design.

Use easily on the go



Express on-the-go without having to worry about a power supply. Our rechargeable battery comes with a micro-USB adapter and lasts up to 3 sessions on a single full charge.

Clean and setup easily



Our closed expression system means milk stays out of tubing so there's less to clean. Fewer parts also make it a breeze to put back together again.

Highlights

Save your favourite settings



Know what you like? The breast pump automatically remembers your last settings, so all you need to do is sit down and press start.

Track your pumping time

Track how long you express on each breast, and your total pumping time with the digital timer display.

Pause whenever you like



If you need to adjust or take a break, the start/pause button is always at your finger tips.

Breast pump bag and pouch



Use the tailor-made bag and insulation pouch to take your pump with you and minimise bulk.

Specifications

Power

Mains voltage: 100 - 240 V

Material

Bottle: BPA free*, Polypropylene

Teat: BPA free*, Silicone

Breast pump: BPA-free* (food contact parts only)

What is included

Travel bag: 1 pcs

Disposable breast pads: 2 pcs

Expression kit: 2 pcs

Insulation pouch: 1 pcs

Micro-USB adapter: 1 pcs

Motor unit (Rechargeable): 1 pcs

Sealing disc: 2 pcs

4oz Bottle: 2 pcs

Slow Nipple: 2 pcs

Ease of use

Breast pump use: Easy cleaning & assembly, quiet motor, memory function, rechargeable battery, motor display, travel bag & pouch

Functions

No leaning forward: Sit in a comfortable position

Settings: 16 expression levels, 8 stimulation levels

Soft & adaptive cushion: Gentle stimulation



* Based on milk flow initiation time (time to Milk Ejection Reflex - MER) results from clinical trial with 20 participants (Netherlands, 2019) compared to time to MER results for other Philips predecessor pump technology from Feasibility study with 9 participants (Netherlands, 2018)

* 1) Mangel et al. Breastfeeding difficulties, breastfeeding duration, maternal body mass index, and breast anatomy: are they related?. Breastfeeding Medicine, 2019, (109 participants, Israel); (2) Ziemer et al. Skin changes and pain in the nipple during the 1st week of lactation.

* Journal of Obstetric, Gynecologic & Neonatal Nursing, 1993, (20 Caucasian participants, USA); (3) Ramsay et al. Anatomy of the lactating human breast redefined with ultrasound imaging, 2005, (28 participants, Australia).

* Based on questionnaire results for 1k cushion from clinical trial with 20 participants (Netherlands, 2019)

* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011