



Philips VisaCare
Microdermabrasion

Skin resurfacing treatment
with Sensitive and Normal tip

SC6240/02

Firmer feeling, more youthful looking skin

with Philips VisaCare Microdermabrasion System

Invigorate and maintain your skin's natural youthful look with VisaCare's Dual Action Air Lift & Exfoliation system. Massages the skin, encouraging blood circulation, while removing dead skin cells, allowing for natural cell renewal.

Microdermabrasion technology

- The Air Lift System gently stimulates microcirculation

Safe and gentle, clinically tested

- Developed with expert scientists and dermatologists

PHILIPS

Specifications

Benefits

- Dead skin cells removal: thanks to the exfoliating tip
- Stimulated microcirculation: lead by Air Lift System™
- Epidermal thickening: driven by abrasive system
- Skincare creams: better absorbed, experienced by 87% of users

Ease of use

- Cordless: up to 6 uses without charging
- Battery indicator: illuminated icon
- Customized tips: Normal and sensitive skin
- Easy handling: Ergonomic and elegant design
- Exclusive and convenient stand: charging and storing stand

Power

- Voltage: 100-240 V

- Charging time: 8 hours
- Power system: Rechargeable battery
- Running time: 6 uses of 5 minutes each

Items included

- Cleaning brush: To clean VisaCare tips
- Instructions for use: Full user manual
- Customized tips: Normal and sensitive skin
- Power adapter: 100 - 240 V adapter
- Storage pouch: For storage and travel
- Exclusive and convenient stand: charging and storing stand

Service


- Warranty: 2-year limited warranty

Application areas

- Face and neck

Highlights

Air Lift System



The vacuum massage intensively stimulates microcirculation. Skin is gently lifted and stretched, massaging its deeper layers. Blood flow increases, bringing oxygen and nutrients to the surface. For a revitalized look. VisaCare is designed to be in the sweetspot of stretch required to stimulate collagen and elastin production. For a firmer feel in 4-6 weeks time.

Developed with skin experts

Philips VisaCare uses microdermabrasion technology commonly used in professional beauty salons. Philips has adapted this technology for safe use in the comfort of your own home. Philips VisaCare has been developed in collaboration with expert scientists and dermatologists.



Issue date 2022-06-16

Version: 3.1.1

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Independent test results on 75 consumers at 4 weeks treatment, UK 2014