

# PHILIPS

## sonicare

2-pack sonic electric toothbrushes - white

3100 series

Pressure sensor

2 minute SmartTimer

Up to 3 x better plaque removal

14 day long-lasting battery life



HX3675/13



Goodbye manual toothbrush. Hello Sonic technology.

Up to 3 x better plaque removal\*

Sonic technology combined with our brush action gently removes plaque up to 3 x better\* than a manual toothbrush.

#### **Designed around you**

- Always know when to replace your brush heads
- 14 day long-lasting battery life
- Ergonomic design makes the toothbrush easy to hold and use
- Goodbye adapter. Hello USB.

#### **Optimise your brushing**

- Optimise your brushing with SmarTimer and QuadPacer

#### **Easy transition to Sonicare**

- Easy-Start to ease the transition

#### **Protects your gums**

- Built-in pressure sensor protects your teeth and gums

#### **Provides superior clean**

- Up to 3 x better plaque removal vs. a manual toothbrush\*
- Our unique technology brings you a powerful yet gentle clean


# Highlights

## Advanced Sonicare Technology




Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two months' worth of manual brushing in just 2 minutes.\*\* 31,000 brush strokes per minute gently clean your teeth, break up plaque and sweep it away for an exceptional daily clean.

## Removes plaque 3 x better\*



The Sonicare electric toothbrush with advanced sonicare technology is clinically proven to remove plaque up to 3 X better\* than a manual toothbrush. It removes plaque from your teeth and along the gum line while protecting your gums.

## Protects your gums



Built-in pressure sensor automatically detects the pressure you apply, warns you and reduces the toothbrush vibrations automatically to help

protect your gums. The toothbrush will make a pulsing sound as a reminder to ease off the pressure. 7 out of 10 people found this feature helped them become better brushers.

## SmartTimer and QuadPacer



2 minute SmartTimer and 30 second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

## Designed around you



Slim, lightweight and ergonomic design makes the toothbrush easy to hold and use so that you get an amazing clean without any strain.

## BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync

technology tracks how long you've been using your brush head for and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep remind you when it's time to replace it.

## Easy-Start



Our Easy-Start programme gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

## Long-lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

## Goodbye adapter. Hello USB.

Power adapters may seem small, but added together, their plastic makes a huge impact on the environment. By switching to a USB charger in our Philips Sonicare 1000-4000 Series alone, we can use up to 189,210 kg less plastic every 4 years. That's powerful, right?

# Specifications

## Power

Voltage: DC5V

## Technical specifications

Operating time (full to empty): 14 days\*\*\*\*\*

Battery: Rechargeable

Battery type: Lithium ION

## Design and finishing

Colour: White

## Service

Warranty: 2-year limited warranty

## Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

BrushSync Replacement Reminder:

Replacement reminder. Always know when to, replace brush heads, Reminder icon lights up, to always ensure best results

## Items included

Brush heads: 2 C1 ProResults brush heads

Charger: 2

Handle: 2 3100 Series

## Cleaning performance

Speed: Up to 31,000 brush strokes per minute

Performance: Removes up to 3 x more plaque\*

Timer: Quadpacer and Smartimer

Pressure feedback: Vibrates handle to alert user

## Modes

Clean: For exceptional everyday cleaning

## Smart sensor technology

BrushSync Replacement Reminder: Always know when to, replace brush heads

Pressure sensor: Alerts when brushing too hard

