

PHILIPS

sonicare

Interdental -
Rechargeable

AirFloss

w/ 2 Nozzles

1 Nozzle holder



HX8233/02



Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is the easiest way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to remove plaque where brushing can't.

Improves oral health

- Removes more plaque between teeth
- Helps improve gum health in just two weeks
- Helps prevent cavities between teeth

Achieve a 30-second clean

- Cleans your entire mouth in just 30 seconds
- Point. Press. Clean!

Easy to use

- Air and micro-droplet technology
- Comfortable ergonomic handle
- Guidance tip ensures correct placement
- A nozzle for each member of the family

Helps build a healthy habit

- An easy way to improve cleaning between teeth


Highlights

An easy way to start




Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

Removes more plaque between teeth



Removes more plaque between teeth

Air and micro-droplet tech



AirFloss utilises unique air and micro-droplet technology to send a burst of pressurised air and water or mouthwash between your teeth. The one-touch activation button is as simple as a mouse click.

Guidance tip



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

30-Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

Ergonomic handle

Designed to be easy to manoeuvre in tight spaces, AirFloss has a comfortable ergonomic handle and a non-slip grip.

Coloured nozzle tips

Coloured nozzle tips and a handy nozzle holder make it easy for the whole family to use AirFloss.

Specifications

Power

Voltage: Multi-voltage charger

Technical specifications

Battery: Rechargeable

Operating time: Replace nozzle every 6 months

Battery type: Lithium ION

Charge time: 24 hours to full charge

Design specifications

Colour(s): White with green accents

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Battery indicator: Flashing indicates low battery

Battery Life: Up to 2 weeks

Nozzle attachment: Easily snaps on and off

Items included

Charger: 1 multi-voltage charger

AirFloss Handle: 1

AirFloss Nozzle: 2

Nozzle holder: 1

Cleaning performance

Auto-burst feature: Cleans entire mouth in 30 secs

Health benefits: Improves gum health in only two weeks

Nozzles: Best results change every 6m

Performance: Removes up to 5x more plaque*

