



Philips Sonicare Power
Flosser 7000 System
Oral Irrigator System

Quad Stream technology

Whisper quiet design
Pulse waves guide as you go
with ExpertClean 7300



HX3921/42

Up to 9X as effective in removing plaque*

An effortless, thorough clean

Enjoy an effortless, thorough clean between your teeth and improve your gum health. Our Quad Stream technology helps you cover more area with less effort, while Pulse Wave technology gently guides you from tooth to tooth. Easy.

An effortless clean with advanced technology

- Easy to fill 600ml reservoir
- Quad Stream technology for faster, more effective flossing
- Standard nozzle for a regular clean
- Pulse Wave technology guides from tooth to tooth
- Whisper quiet for great results without disturbing others
- Quad Stream nozzle for maximum coverage
- Compact nozzle case
- Philips Sonicare ExpertClean 7300
- Easy magnetic dock with auto shut-off

Customize easily to your needs


- 4 flossing modes, 10 intensities

PHILIPS

sonicare


Highlights

Covers more area, effortlessly



X-shaped water streams cover more surface area for a faster, deeper clean versus manual floss. Effortlessly.

Guides from tooth to tooth



Gentle pulses of water guide you from tooth-to-tooth so that you don't miss a spot.

Whisper quiet design



Engineered for a discreet clean, so you get great results without disturbing others.

4 modes, 10 intensities




Enjoy a clean that's tailored to your needs. Clean mode uses a continuous flow of water for a great everyday clean; Deep Clean uses distinct pulses for a more thorough clean; Sensitive mode is extra gentle for sensitive gums; Gum care mode soothes and stimulates the gumline. Whenever you need, adjust the intensity for your comfort using the LED light display as your guide.

F3 Quad Stream nozzle



For maximum coverage and an extra thorough clean, use the Quad Stream nozzle to concentrate water into four wide streams between teeth and along the gum line. The soft rubber guide is gentle on gums and every nozzle simply clicks on and off.

F1 Standard nozzle



Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth.

Magnetic dock w/auto shut-off



When you're finished, simply return the handle to the easy magnetic dock. The handle automatically switches from pause to off.

600ml reservoir



The 600ml reservoir with easy flip lid holds enough water for one cleaning. To refill simply lift out, fill, and click back in.

ExpertClean 7300



Complete your clean with the Philips Sonicare ExpertClean 7300. It's designed to give you an exceptional daily clean and help you keep up good brushing habits. Smart sensors help you maintain the perfect pressure and an in-app progress report encourages you to keep up your routine. When you're done, simply return the handle to the built-in charger, or use the travel charger when you're on the go.

Compact nozzle case

Store extra nozzles hygienically in our compact and stylish nozzle case.

Specifications

Power

- Voltage: Multi-voltage charger

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

Performance

- Cleaning: Clean entire mouth in 60-90sec

Items included

- Brush heads: 1 C3 Premium Plaque Control
- Charger: 1
- Travel case: 1

- Handle: 1 ExpertClean
- F1 Standard nozzle: 1
- F3 Quad Stream nozzle: 1
- Power Flosser: 1
- Nozzle storage case: 1

Modes

- Clean: For exceptional everyday clean
- Gum Care: Gently massages gums
- Deep Clean+: For an invigorating deep clean
- Intensities: 10
- Sensitive: Gentle teeth and gum cleaning

Health benefits

- Plaque removal: 9x as effective*
- Gum health: Up to 270% more effective**



Issue date 2023-09-17

Version: 3.3.1

EAN: 87 10103 94729 5

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* in hard-to-reach areas vs. manual toothbrush + floss. Using Sonicare Power Flosser with Quad Stream tip on setting 8 and Sonicare ExpertClean power toothbrush with Premium Plaque Control brush head in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss

** vs. manual toothbrush + floss. Using Sonicare Power Flosser with Quad Stream tip on setting 8 and Sonicare ExpertClean power toothbrush with Premium Plaque Control brush head in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss