



Philips Sonicare Cordless
Power Flosser 2000
Oral Irrigator

Removes up to 99% of plaque in treated areas*

A thorough clean in-between

Effortlessly thorough cleaning between teeth to improve your gum health. Adjustable intensity and 360° rotation adapts to your needs to clean away plaque and food debris in hard-to-reach areas.

Easy and convenient

- Easy-to-fill 250 ml reservoir
- Fast, universal charging with USB-C cable
- Cordless, compact design makes flossing simple
- Powerful and gentle for a thorough clean

Adapts to your needs

- 3 intensity settings to suit your needs
- 360° rotation for hard-to-reach spots



HX3801/11

PHILIPS
sonicare

Specifications

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

Performance

- Cleaning: Clean entire mouth in 60-90sec

Items included

- Cordless Power Flosser: 1

- F1 Standard nozzle: 1
- USB charging cable: 1

Modes

- Clean: For exceptional everyday clean
- Intensities: 3

Health benefits

- Plaque removal: Removes up to 99% of plaque in treated areas*

Highlights

A powerful, gentle clean

Designed to remove plaque and food debris even from hard-to-reach areas. Leave your mouth feeling fresh and clean without the hassle of string flossing.

Cordless, compact design

Compact, ergonomic design offers a convenient solution for everyday flossing. Cordless handle allows for easy portability and flossing anywhere.

360° nozzle rotation



Access those hard-to-reach areas with a nozzle that rotates 360 degrees, so that you can use the device in any orientation.

3 intensity settings


Choose the clean that suits your smile. Adjustable intensity settings adapt for comfort and cleaning power.

250ml reservoir



Reservoir capacity of 250 ml holds enough water for a recommended 60-second clean. Simply twist and remove or use the side-fill door for easy refilling.

USB-C Charging



Convenient charging with a USB-C cable uses the same connection as most mobile phone or tablet. Get up to 14 days of cleaning from a single charge.



Issue date 2023-09-01

Version: 3.3.2

EAN: 87 10103 99294 3

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* In an in-vitro study. Actual results may vary.