

Philips Sonicare Cordless  
Power Flosser 2000  
Oral Irrigator



HX3801/11

Removes up to 99% of plaque in treated areas\*

A thorough clean in-between

Effortlessly thorough cleaning between teeth to improve your gum health. Adjustable intensity and 360° rotation adapts to your needs to clean away plaque and food debris in hard-to-reach areas.

**Easy and convenient**

- Easy-to-fill 250 ml reservoir
- Fast, universal charging with USB-C cable
- Cordless, compact design makes flossing simple
- Powerful and gentle for a thorough clean

**Adapts to your needs**

- 3 intensity settings to suit your needs
- 360° rotation for hard-to-reach spots

**PHILIPS**  
**sonicare**

# Highlights

## A powerful, gentle clean



Designed to remove plaque and food debris even from hard-to-reach areas. Leave your mouth feeling fresh and clean without the hassle of string flossing.

## Cordless, compact design



Compact, ergonomic design offers a convenient solution for everyday flossing. Cordless handle allows for easy portability and flossing anywhere.

## 360° nozzle rotation



Access those hard-to-reach areas with a nozzle that rotates 360 degrees, so that you can use the device in any orientation.

## 3 intensity settings



Choose the clean that suits your smile. Adjustable intensity settings adapt for comfort and cleaning power.

## 250ml reservoir



Reservoir capacity of 250 ml holds enough water for a recommended 60-second clean. Simply twist and remove or use the side-fill door for easy refilling.

## USB-C Charging



Convenient charging with a USB-C cable uses the same connection as most mobile phone or tablet. Get up to 14 days of cleaning from a single charge.

# Specifications

## Consumer Trade Item

- Height: 24.7 cm
- Width: 13.4 cm
- Length: 9.2 cm
- Net weight: 0.38 kg
- Gross weight: 0.598 kg
- GTIN: 8710103992943
- Country of origin: CN
- Harmonized system code: 85098000

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

## Performance

- Cleaning: Clean entire mouth in 60-90sec

## Items included

- Cordless Power Flosser: 1
- F1 Standard nozzle: 1
- USB charging cable: 1

## Modes

- Clean: For exceptional everyday clean
- Intensities: 3

## Health benefits

- Plaque removal: Removes up to 99% of plaque in treated areas\*



Issue date 2024-08-27

Version: 3.3.2

EAN: 87 10103 99294 3

© 2024 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* In an in-vitro study. Actual results may vary.