

PHILIPS

sonicare

Sonic electric
toothbrush - Dispense

EasyClean

1 mode

1 brush head

HX6581/50



Advanced plaque removal made simple

The switch from a manual toothbrush made easy

Philips Sonicare EasyClean combines our patented Sonicare toothbrush technology and the ProResults brush head to remove up to 2x more plaque than a manual toothbrush. One-mode simplicity and a slim handle for comfort and maneuverability.

Designed around you

- EasyStart
- Recharge indicator

Guides you to follow dental recommendations

- Quadpacer
- Smartimer

Provides a superior clean

- Optimized Philips Sonicare performance
- Patented Sonicare toothbrush technology


Highlights

Sonic technology



Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

EasyStart




Gently increases power over the first 14 uses to ease patients into the Sonicare brushing experience.

Smartimer



Two-minute timer helps ensure dental professional recommended brushing time.

Quadpacer



Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

Optimized performance



This Philips Sonicare toothbrush head features a contoured profile to naturally fit the shape of your teeth and to clean hard-to-reach areas.

Recharge indicator

Light blinks to let user know the battery needs recharging.

Specifications

Power

Voltage: 110-220 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 2 weeks

Battery type: Lithium ION

Design and finishing

Color: Glacier Green

Service

Warranty: 2-year limited warranty

Ease of use

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design

Brushing time: Up to 2 weeks

Items included

Handles: 1 EasyClean

Charger: 1

Brush heads: 1 ProResults standard

Cleaning performance

Performance: Up to 2X more plaque removal*

Health benefits: Helps improve gum health

Speed: Up to 62000 brush movement/min

Timer: Quadpacer and SmarTimer

Whitening benefits: Whitens teeth up to 2 shades

Modes

Clean: For exceptional everyday clean

