

# PHILIPS

## sonicare

Sonic electric  
toothbrush

3 Series gum health

1 mode

3 intensity settings

2 brush heads



HX6632/18



## Improves gum health up to 100%\*

\*more than manual brushing

Provides superior plaque removal and gum care in an easy-to-use power toothbrush to help prevent and reduce the early signs of gum disease while offering your most comfortable brushing experience

### Guides you to follow dentist recommendations

- SmarTimer helps you meet recommended brushing time
- Quadpacer interval timer encourages thorough brushing

### Improves gum health better vs. a manual toothbrush

- Better than a manual toothbrush for gingivitis
- Removes up to 6X more plaque along the gumline
- Sonicare dynamic cleaning action drives fluid between teeth

### Offers a customized brushing experience

- Three intensity settings for a gentle experience

### Recommended by dental professionals world wide

- Proven to improve gum health 3 ways
- Safe on orthodontics, dental work and veneers

# Highlights

## Improves gum health 3 ways




Clinically proven to improve gum health three ways: reduces redness, swelling, bleeding (sites)

## Prevents & reduces gingivitis



With superior plaque removal along the gum line, this toothbrush can help reduce gingivitis up to 100% better than a manual toothbrush in just two weeks

## Three intensity settings



With three customizable intensity settings (low, medium & high) you can select the right power level to ensure that brushing along the gumline is comfortable. After 1.5 seconds of brushing, you will have to push twice on button to switch mode. Once to pause and second to switch.

## Up to 6x better plaque removal



Removes up to 6X more plaque along the gumline than a manual toothbrush.

## Sonic technology



Our sonic technology creates a unique dynamic cleaning action that drives fluid deep between teeth and along the gumline for a thorough-yet-gentle clean

## Safe and gentle




The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worry.

## Smarter



This Philips Sonicare power toothbrush features a special Smarter to help you meet the dentist-recommended brushing time of two minutes

## Quadpacer



Features a Quadpacer interval timer to help encourage thorough brushing of the entire mouth. The toothbrush beeps every 30 seconds, signifying it's time to move to another area of the mouth.

# Specifications

## Power

Voltage: 110-220 V

## Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 3 weeks\*\*

Battery type: Lithium ION

## Design and finishing

Color: Persian Green

## Service

Warranty: 2-year limited warranty

## Ease of use

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design

Brushing time: Up to 3 weeks\*\*

## Items included

Handles: 1 3 Series gum health

Charger: 1

Travel case: 1

Brush heads: 2 ProResults gum health

Mouth rinse: 27 ml Listerine

## Cleaning performance

Performance: Removes up to 6X more plaque\*

Health benefits: Improves gum health

Speed: Up to 31000 brush movement/min

Timer: Quadpacer and SmarTimer

## Modes

Clean: For exceptional everyday clean

3 intensities: Low, Medium, High

