

PHILIPS

sonicare

Sonic electric
toothbrush

EasyClean

1 mode

2 brush heads

BreathRX and tongue scraper

Travel bag

HX6511/34



Better plaque removal

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

Proven to improve oral health

- Clinically proven safe and gentle
- Philips Sonicare toothbrush helps whiten teeth

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- Removes up to 2x more plaque than a manual toothbrush
- Dynamic cleaning action drives fluid between teeth

Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience

Highlights

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

SmarterTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Helps whiten teeth



This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

Patented Sonic technology



Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.



Specifications

Design and finishing

Color: White

Power

Voltage: 110-220 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 2 weeks

Battery type: Lithium ION

Service

Warranty: 2-year limited warranty

Ease of use

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design

Brushing time: Up to 2 weeks

Items included

Handles: 1 EasyClean

Charger: 1

Brush heads: 1 gum health brush head, 1 ProResults standard

Travel bag: 1

Tongue scraper: 1

BreathRX mouth rinse (2oz.): 1

Cleaning performance

Performance: Up to 2X more plaque removal*

Health benefits: Helps improve gum health

Speed: Up to 62000 brush movement/min

Timer: Quadpacer and SmarTimer

Whitening benefits: Whitens teeth up to 2 shades

Modes

Clean: For exceptional everyday clean

