

# PHILIPS

## sonicare

Sonic electric  
toothbrush

For Kids

2 modes

2 brush heads



HX6311/02



## Encourage healthy brushing on their own Electric toothbrush for kids

The Philips Sonicare electric toothbrush HX6311/02 supports growing smiles

### **Provides a superior clean**

- Sonicare dynamic cleaning action drives fluid between teeth

### **Clean and protect growing smiles**

- Rubberized brush head is designed to protect young teeth

### **Designed to develop better brushing habits**

- Age-appropriate brush heads to protect kid's teeth
- Anti-roll shape
- Multi grip design for parents and kids
- 2 Kid-friendly power modes deliver a gentle, effective clean
- Fun interchangeable panels for customization
- KidPacer delivers fun and engaging musical tones
- KidTimer helps to increase brushing time

# Highlights

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Rubberized brush head



Rubberized brush head is designed to protect young teeth

## KidTimer




Increases brushing time slowly over 90 days, building towards the dentist recommended 2 minutes

## Power modes




2 Kid-friendly power modes adapts to different ages providing a gentle cleaning that's ideal for different ages. There is a low mode for ages 4 and up and high mode for ages 7 and older

## Age-appropriate brush heads



This Philips Sonicare electric toothbrush has 2 brush head sizes that are specifically designed to gently clean and protect teeth at key developmental stages

## KidPacer



KidPacer encourage kids to brush their entire mouth effectively for the dentist-recommended 2 minutes

## Interchangeable panels



Allows kids to customize their toothbrush

## Anti-roll shape



Sturdy handle design allows kids to store upright as well as apply toothpaste while the toothbrush is laying flat.

## Multi grip design



Ergonomically designed to enable parents and kids to brush together or own their own

# Specifications

## Power

Voltage: 110-220 V

## Technical specifications

Operating time (full to empty): Up to 3 weeks\*\*

Battery: Rechargeable

Battery type: Lithium ION

## Design and finishing

Color: Red with 3 interchangeable panels

## Service

Warranty: 2-year limited warranty

## Ease of use

Display: Illuminated display

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design, Rubber grip for easy handling

Brushing time: Up to 3 weeks\*\*

## Items included

Charger: 1

Handles: 1 Sonicare for Kids

Brush heads: 1 Sonicare for Kids standard, 1 Sonicare for Kids compact

## Cleaning performance

Performance: 75% more effective\*

Health benefits: For healthy oral care habits

Speed: Up to 62000 brush movement/min

Timer: KidTimer and Quadpacer

## Modes

Power modes: 2

