

PHILIPS

AVENT

DECT Baby Monitor

Audio Monitors

100% private connection

Night light and lullabies

Talkback function

Vibration alert

 Product Picture

SCD570/01



The most reliable connection to your baby

Advanced monitoring for you and your baby

Our new DECT SCD570/01 Baby Monitor delivers advanced monitoring. Providing the most reliable connection with crystal clear sound, temperature sensor, a soothing night light and lullabies for you and your baby.

Freedom and flexibility

- Range up to 330 metres*
- A docking station to recharge your parent unit
- Superior operating time up to 18 hours
- Vibration alert
- Instant dim mode for a peaceful night

Comforting and soothing your baby

- Stay connected with talk back
- Comforting night light and lullabies to calm your little one

Reliable, secure connection and constant reassurance

- Always know when your monitor is in range and connected
- Perfect clear sound thanks to DECT Technology
- Energy-saving Smart ECO mode
- Sound-level lights indicate if the baby makes a noise
- Monitor the temperature in your baby's room

Highlights

Zero interference with DECT



DECT technology guarantees zero interference from any other transmitting products like other baby monitors, cordless phones and mobile phones. Data encryption provides a secure and private connection, so you can be sure that you are the only one that can hear your baby.

Perfect clear sound



Hear every tiny giggle, gurgle and hiccup with perfect clarity. DECT (Digital Enhanced Cordless Telecommunications) technology provides high quality, crystal clear sound so you can hear your baby all the time.

Smart ECO mode



The unique Smart ECO mode automatically minimises your transmission power and increases your battery lifetime. The closer you


are to your baby, the less power is needed for a perfect connection (not available in US and Canada).

Sound activation 5 LED lights




Lights are activated to indicate the level of sound in your baby's room, even when the Parent Unit is muted.

Always connected



The parent unit will alert you when your monitor is out of range or when the power is low, helping you to make sure that you are always connected to your baby.

Temperature sensor



The temperature sensor with customised alert allows you to monitor the climate in your baby's room. Your baby can't regulate its body temperature as well as you can and a slight change in temperature can make your baby restless. A customised temperature sensor will

immediately alert you, through the digital display on the parent unit, if the climate in your baby's room changes.

330 m Range



Indoor range up to 50 metres*. Outdoor range up to 330 metres*.

Docking station

A docking station enables you to conveniently recharge your parent unit.

Operating time up to 18 hours



The small, rechargeable parent unit will give you the freedom of cordless streaming for up to 18 hours before re-charging is required.

Vibration alert

Personalise your alert settings and be notified by a silent and subtle vibration of the parent unit.

Peaceful night mode

Press the night dim button to instantly activate the peaceful night mode. In this mode the light intensity of the lights, the display, as well as the volume of the alerts are reduced. Easily press the night dim button again in order to switch back to normal operating mode.

Specifications

Tuner/Reception/Transmission

Frequency band: DECT
Automatic channel selection

Convenience

Auto out-of-range warning
Battery low indication
Battery charging indication
Volume control
Sensitivity control
Sound-level lights

Power

Operating time on battery hour: 18 hours
Charging time: 10 hour(s)
Power Supply: 120 V (US), 220–240 V

Technical specifications

Storage temperature range: 10–40 °C
Operating temperature range °C: 10 – 40 °C

Features

Comforting night light
Charging station
DECT Technology
Indoor range up to: 50 m
Energy-saving Smart ECO mode: ECO mode,
Yes

Linked and in range indication

Outdoor range up to: 330 metres
Talkback function
Temperature sensor
Peaceful night mode and vibration alert
LCD Screen
Sound activation lights: 5 LED lights
Soothing lullabies

Logistic data

F-box dimensions (W x H x D): 220 x 156 x
94 mm

Development stages

Stage: Pregnancy, 0–6 months

