

PHILIPS

sonicare

Rechargeable
toothbrush

1600-Series

1 mode

HX1630/02



Gentle on gums, tough on plaque

The ultimate electric toothbrush

Gums are an important part of your oral health. That's why the new Philips rechargeable electric toothbrush HX1630 features a patented Gum Protection system that adjusts to maintain an optimal brushing pressure.

Double Cleaning Action

- Naturally whiter teeth
- Targets hard-to-reach plaque in between teeth
- Cleans visible teeth surfaces

Gentle on gums, tough on plaque

- Gum Protection system
- Adjusts to maintain optimal brushing pressure

Additional model features

- 2-minute timer
- 1 extra free brush head! 2 in total
- Recharge indicator

Highlights

Gum Protection system



Philips Sonicare's patented Gum Protection system adjusts to maintain optimal brushing pressure

Patented Gum Protection system



Adjusts to maintain optimal brushing pressure

Double Cleaning Action



Cleans visible teeth surfaces while the active tip targets hard-to-reach plaque in between teeth.

Targets

Targets hard-to-reach plaque in between teeth

Naturally whiter teeth



Naturally whiter teeth

1 extra free brush head!


1 extra free brush head! 2 in total

2-minute timer



Dentists' recommend a 2 minute brushing time.

Recharge indicator



Lights up when charging properly.

Specifications

Technical specifications

Cord length: 1.16 m

Power consumption: 2.9 W

Voltage: 220-240 V

Frequency: 50/60 Hz

Charging time: 16 hour(s)

Operating time (full to empty): 25 minute(s)

Dual brush head: Brush head combines to create 15,000 strokes/min.

Battery rechargeable: NiCd

Weight and dimensions

A-box dimensions: 261 x 247 x 247 mm

F-box weight: 339

A-box number of pieces: 6

F-box dimensions (L x W x H): 117 x 82 x 225 mm

A-box weight: 2354 g

Accessories included

Charger base with brush head storage pins

Brush head: 1

Additional features

1 extra free brush head

Two minute timer

Recharge indicator

