

PHILIPS

sonicare



Trial Unit

ProtectiveClean 4300

Built-in pressure sensor

1 cleaning mode

1 x BrushSync feature



HX6848/98

Healthier gums. Gently does it.

Improves gum health up to 100% more than a manual

Dental Professional Trial Unit for dental professional staff to try the power of Sonicare, with exclusive discounted pricing on the ProtectiveClean 4300

Proven to improve oral health

- Philips Sonicare's advanced sonic technology

Innovative technology

- A gentle nudge to go nice and easy
- Patients get the most from their brush heads
- Connects smart brush handle and smart brush heads

Improves patients' oral health

- Improves gum health up to 100%*
- Fillings, crowns, or sensitive areas? Not to worry.

Designed around your patient

- A great way to get used to your new toothbrush
- Encouragement to brush thoroughly

Highlights

Improves gum health



Densely packed, high-quality bristles give you an extra gentle brushing to remove plaque along the gum line to improve gum health up to 100% more than a manual toothbrush. And the specially curved power tip makes reaching the teeth at the back of your mouth a breeze.

Too much pressure?



Patients may not notice if they are brushing too hard, but their toothbrush will. If patients start taking it out on their gums, the toothbrush will make a pulsing sound as a reminder to ease off on the pressure.

Safe and gentle




You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns and veneers, and helps prevent cavities and improve gum health.

Philips Sonicare technology




Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

Smart replacement reminders



All brush heads wear out over time. But we've got that covered too: our smart brush-recognition technology tracks how long patients have been using their brush head for and how hard they've been brushing. When it's time to replace it, a light on the handle lets them know. That way, they can be sure that their brush head is consistently doing a stellar job.

BrushSync technology




A microchip-enabled technology that detects and synchronises the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

Helpful timers



Our QuadPacer lets patients know when they've spent just the right amount of time cleaning each part of their mouth, while our Smartimer tells them when they've brushed for the recommended two minutes.

Easy-start programme



Starting something new can take getting used to. That's why our easy-start programme gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

ProtectiveClean 4300

Trial Unit

HX6848/98

Specifications

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Colour: Mid-Blue

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

Replacement reminder: To always ensure best results, reminder icon lights up

Items included

Brush heads: 1 G3 Premium Gum Care, 1 InterCare standard, 1 W2c Optimal White compact

Charger: 1

Handle: 1 ProtectiveClean

Cleaning performance

Health benefits: Improves gum health in only two weeks, Helps reduce cavities

Speed: Up to 62,000 brush movements/min

Timer: Quadpacer and Smartimer

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday cleaning

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

