

# PHILIPS

## sonicare

Sonic electric  
toothbrush

ProtectiveClean  
5100

Built-in pressure sensor

3 modes

1 x BrushSync feature

Travel case



HX6856/10



## Whiter teeth. Gently does it.

Whitens teeth in just 1 week.

Feel the difference of a gentle clean with our Pressure Sensor while whitening your teeth in 1 week.

### Proven to improve oral health

- Improves gum health by up to 100%\*
- Safe and gentle on sensitive areas, orthodontics and dental work
- Whitens teeth in just one week

### Designed around you

- Encouragement to brush thoroughly
- Makes traveling easier for you
- Always know when to replace your brush heads

### Innovative technology

- Lets you know if you're pressing too hard
- Philips Sonicare's advanced sonic technology
- Connects smart brush handle and smart brush heads

### Optimise your brushing

- Choose from three modes


# Highlights

## Naturally whiter teeth



Click on the W2 Optimal White brush head to remove surface stains and reveal a whiter smile. With its densely packed central stain-removal bristles, it's clinically proven to whiten teeth in just one week.

## Improves gum health



Densely packed, high-quality bristles give you an extra gentle brushing to remove plaque along the gum line to improve gum health up to 100% more than a manual toothbrush. And the specially curved power tip makes reaching the teeth at the back of your mouth a breeze.

## Three modes



With a choice of three modes, this toothbrush allows you to customise brushing to your own needs. Clean mode is standard for superior cleaning. White is the ideal mode for removing surface stains. And Gum Care mode

adds an extra minute of reduced-power brushing, so you can gently massage your gums.

## Safe and gentle



You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns and veneers, and helps prevent cavities and improve gum health.

## Monitor your brushing pressure



Brushing too hard can damage your teeth and gums. To prevent this, your Philips Sonicare ProtectiveClean emits a gentle pulsing sound to remind you to ease off on the pressure.

## BrushSync technology




A microchip-enabled technology that detects and synchronises the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

## BrushSync replacement reminder



All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

## Helpful timers



Need an electric toothbrush with a timer? Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

## Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

# Specifications

## Design and finishing

**Colour:** Pastel pink

## Power

**Voltage:** 110-220 V

## Technical specifications

**Operating time (full to empty):** Up to 2 weeks

**Battery:** Rechargeable

**Battery type:** Lithium ION

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Handle:** Slim ergonomic design

**Handle compatibility:** Easy click-on brush heads

**Battery indicator:** Light shows battery status

**Brush head recognition:** Syncs with optimal mode

**Brushing time:** Up to 2 weeks

**Replacement reminder:** To always ensure best results, reminder icon lights up

## Items included

**Brush heads:** 1 W2 optimal white, 1 G2 Optimal Gum Care

**Travel case:** 1

**Handle:** 1 ProtectiveClean

**Charger:** 1

## Cleaning performance

**Speed:** Up to 62,000 brush movements/min

**Health benefits:** Helps reduce cavities, Improves gum health

**Whitening benefits:** Whitens teeth in 1 week

**Timer:** Quadpacer and Smartimer

**Pressure feedback:** Vibrates handle to alert user

## Modes

**Clean:** For exceptional everyday cleaning

**Gum Care:** Gently massages gums

**White:** Removes surface stains

## Smart sensor technology

**BrushSync Replacement Reminder:** Always know when to, replace brush heads

**Pressure sensor:** Alerts when brushing too hard

