

# PHILIPS

## sonicare

Sonic electric  
toothbrush

FlexCare+

5 modes

2 brush heads

HX6922/03



## Superior cleaning between teeth

Toothbrush for gentle cleaning of teeth and gums

Combining our sonic technology and 5 brushing modes, FlexCare+ electric toothbrush is your solution for improved oral health. Plus you can improve your gum health after just 2 weeks of use.

### Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology
- Improves gum health in only two weeks

### Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

### Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- Designed to reach deeper and go further

### Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience
- Five different modes for a better cleaning experience

# Highlights

## Reach deeper and go further



Philips Sonicare InterCare toothbrush head features extra-long, high-density bristles to target hidden plaque caught deep between teeth and other hard-to-clean areas.

## Naturally whiter teeth




The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning along the gumline for improved gum health in just two weeks. Up to 6x better plaque removal than a manual toothbrush at the gumline.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## SmarterTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Five different cleaning modes




The five modes include a Clean mode for an outstanding clean in two minutes. A Sensitive mode for gentle yet thorough cleaning of the gums. A Refresh mode for those quick one minute touch-ups. A Massage mode to stimulate the gum and a Gum Care mode that helps improve gum health.

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth



# Specifications

## Power

Voltage: 110–220 V

## Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 3 weeks\*\*

Battery type: Lithium ION

## Design and finishing

Color: Cooper frost

## Service

Warranty: 2-year limited warranty

## Ease of use

Brush head system: Easy click-on brush heads

Display: Illuminated display

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

Brushing time: Up to 3 weeks\*\*

## Items included

Charger: 1

Brush heads: 2 InterCare standard

Travel case: 1

Handles: 1 FlexCare+

## Cleaning performance

Performance: Removes up to 6X more plaque\*

Speed: Up to 62000 brush movement/min

Health benefits: Helps improve gum health

Whitening benefits: Helps naturally remove stains

Timer: Quadpacer and SmarTimer

## Modes

Clean: For exceptional everyday clean

Gum Care: Gently massages gums

Massage: Invigorating massage

Sensitive: Gentle teeth and gum cleaning

Refresh: 1 min mode for quick touch up

