

# PHILIPS

## sonicare



2-pack sonic electric  
toothbrushes with cases

ProtectiveClean 4300

Built-in pressure sensor

1 cleaning mode

1 x BrushSync feature

Travel case



HX6800/35

## Cleaner teeth. Gently does it.

Remove up to 7 x more plaque than a manual

Feel the difference of a gentle clean with our Pressure Sensor while removing up to 7 x more plaque than a manual toothbrush

### Proven to improve oral health

- Remove up to 7x more plaque with our InterCare brush head\*
- Safe and gentle on sensitive areas, orthodontics and dental work

### Designed around you

- Easy-start programme builds your Philips Sonicare routine
- Encouragement to brush thoroughly
- Makes traveling easier for you

### Innovative technology

- Lets you know if you're pressing too hard
- Always know when to replace your brush heads
- Connects smart brush handle and smart brush heads
- Philips Sonicare's advanced sonic technology

# ProtectiveClean 4300

2-pack sonic electric toothbrushes with cases

HX6800/35

## Highlights

### Remove up to 7 x more plaque



Click on our InterCare brush head to improve gum health in just 2 weeks. Extra-long bristles help remove more plaque in hard-to-reach spots and between teeth for healthy gums.

### Monitor your brushing pressure



Brushing too hard can damage your teeth and gums. To prevent this, your Philips Sonicare ProtectiveClean emits a gentle pulsing sound to remind you to ease off on the pressure.

### Safe and gentle




You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns and veneers, and helps prevent cavities and improve gum health.

### Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

### BrushSync replacement reminder



All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

### Helpful timers




Need an electric toothbrush with a timer? Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

### BrushSync technology



A microchip-enabled technology that detects and synchronises the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

### Ease into your new routine



Starting something new can take some getting used to. Our easy-start programme gives you the option of a gradual, gentle increase in brushing power with the first 14 uses of your new toothbrush.

### Handy travel case and charger



Our premium travel case lets you store your toothbrush hygienically, while our compact charging base keeps you topped up when you're on the go. You do get two weeks of regular use from a single full charge, but the charger's just the ticket for longer trips.

# ProtectiveClean 4300

2-pack sonic electric toothbrushes with cases

HX6800/35

## Specifications

### Power

Voltage: 110-220 V

### Technical specifications

Operating time (full to empty): 14 days\*\*

Battery: Rechargeable

Battery type: Lithium ION

### Design and finishing

Colour: Black, Pastel pink

### Service

Warranty: 2-year limited warranty

### Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Battery indicator: Light shows battery status

### Items included

Brush heads: 2 i InterCare

Travel case: 2

Charger: 1

Handle: 2 ProtectiveClean

### Cleaning performance

Health benefits: Improves gum health in only two weeks, Helps reduce cavities

Speed: Up to 62,000 brush movements/min

Performance: Removes up to 7x more plaque\*

Timer: Quadpacer and Smartimer

Pressure feedback: Vibrates handle to alert user

### Modes

Clean: For exceptional everyday cleaning

2 intensities: Low, High

### Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

