

# PHILIPS

## sonicare

Sonic electric  
toothbrush

DiamondClean

5 modes

2 brush heads

Glass charger

USB travel charger



HX9362/67



## Whiter, healthier teeth for life

### Best Philips Sonicare whitening sonic toothbrush

Sonicare's best whitening in our most elegant Philips Sonicare electric toothbrush. Switch to Sonicare.

#### Proven to improve oral health

- A whiter smile in 1 week with our DiamondClean brush head\*
- Improves gum health in just 2 weeks\*
- Up to 7x more plaque removal than a manual toothbrush
- Safe and gentle for your oral health needs
- Choose from 5 modes including Deep Clean and Sensitive

#### Designed around you

- Timers encourage thorough 2 minutes of brushing
- USB charging travel case and charging glass

#### Provides a superior clean

- Angled neck to help tackle cleaning in tough areas
- Our unique technology brings you a powerful yet gentle clean

# Highlights

## Up to 2x whiter teeth




Click on our DiamondClean brush head to gently yet effectively remove surface stains. Densely packed central stain removal bristles work hard to make your smile 2x whiter in just 7 days.\*

## Say goodbye to plaque




Densely-packed, high-quality bristles remove up to 7x more plaque than a manual toothbrush.

## Look after those gums



With optimal cleaning from your DiamondClean, gums will become healthier in 2 weeks\*. Removing up to 7x more plaque along the gum line than a manual toothbrush\*, you'll get your healthiest smile.

## 5 modes to clean your way



With DiamondClean, you get a truly refreshing clean every day. Our 5 modes cover all your brushing needs; Clean mode – for exceptional daily cleaning, Gum Care – to gently massage gums, Deep Clean – for an invigorating deep clean, Sensitive – for gentle yet effective cleaning of sensitive gums, and White – the ideal mode to remove surface stains.

## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Reach those tricky spots



The unique shape of your DiamondClean toothbrush handle combined with our brush heads means that even tricky spots like back teeth get a thorough clean.

## Philips Sonicare technology



Powerful sonic vibrations whip up your toothpaste into plaque fighting bubbles, and drive them deep between your teeth and along your gum line. At the same time your teeth will experience 62,000 gentle but effective brush movements. You'll get a month's worth of manual brushing in just 2 minutes.

## Design with storage in mind



Your deluxe USB travel case doubles as a charger so you stay topped up on the go. Place your toothbrush in the case and connect to your laptop or plug into a wall socket. Also includes a brush head holder for extra hygiene when traveling. As for charging at home, our sleek charger glass fits stylishly into your bathroom and doubles as a glass for rinsing after brushing. Simply place your brush into the glass. Enjoy two weeks of regular use from a single full charge.

# Specifications

## Design and finishing

**Color:** Pink

## Power

**Voltage:** 100-240 V, 50-60 Hz

## Technical specifications

**Operating time (full to empty):** Up to 3 weeks\*\*

**Battery:** Rechargeable

**Battery type:** Lithium ION

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Brush head system:** Easy click-on brush heads

**Battery indicator:** Illuminated icon indicates battery life

**Handle:** Slim ergonomic design

**Display:** Illuminated display

**Brushing time:** Up to 3 weeks\*\*

## Items included

**Travel case:** USB travel charger

**Handle:** 1 DiamondClean

**Brush heads:** 1 DiamondClean compact, 1 DiamondClean standard

**Glass charger:** 1

## Cleaning performance

**Performance:** Removes up to 7x more plaque\*

**Speed:** Up to 62000 brush movement/min

**Health benefits:** Improves gum health in only two weeks

**Whitening benefits:** Whitens teeth 2X better\*

**Timer:** Quadpacer and SmarTimer

## Modes

**Clean:** For exceptional everyday clean

**White:** Removes surface stains

**Gum Care:** Gently massages gums

**Polish:** Brighten and polish your teeth

**Sensitive:** Gentle teeth and gum cleaning

