

PHILIPS

Sonic electric
toothbrush with app

DiamondClean 9000

Connected brushing made easy

Built-in pressure sensor

Smart brush head recognition

4 modes, 3 intensities

HX9911/76



Whiter, healthier teeth for life

Whiter teeth in just 1 day****

Philips Sonicare's best whitening in our most elegant sonic electric toothbrush. Switch to Philips Sonicare.

Enjoy a stylish and user-friendly design

- Simple to charge and style to boot

DiamondClean users experience improved oral health

- Deepest clean while protecting your sensitive gums
- Up to 7x healthier gums* in just 2 weeks

Let the DiamondClean 9000 adapt to your needs

- Personalized brushing experience
- Lets you know if you're pressing too hard
- BrushSync automatically selects the best mode for you
- Always know when to replace your brush heads

Stay on track with a personalized progress report

- Start and keep up healthy habits


DiamondClean 9000

Sonic electric toothbrush with app

HX9911/76

Highlights

Up to 7x healthier gums*




Focus on improving your gum health with the G3 Premium Gum Care brush head. Its smaller size and targeted gum line bristles make it the perfect way to keep this area gum disease free. People who use this brush head have up to 100% less gum inflammation* and up to 7x healthier gums in just two weeks.*

Premium Gum Care brush head



The G3 Premium Gum Care brush head is designed to help you improve your gum health. Its smaller size and targeted gumline bristles gently and effectively clean the gum line, where gum disease commonly starts out.

Progress Report



The DiamondClean 9000 gives you the guidance you need to improve and maintain healthy brushing habits in between dentist checkups. Built-in smart sensors let you know when you're using too much pressure, and by connecting your

brushing experience to the Sonicare app, a personalized Progress Report helps you stay on track, to see how much you've improved over time.

4 modes, 3 intensities



DiamondClean 9000 comes with Clean, White+, Gum Health and Deep Clean+ to take care of your brushing needs. Clean mode is for exceptional everyday cleaning, White+ is for removing stains, Gum Health provides a gentle yet effective clean for gums and Deep Clean+ gives you an invigorating deep clean. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more-sensitive mouths.

BrushSync mode pairing



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. For an example, say you're using the W3 Premium White brush head. Your DiamondClean 9000's BrushSync technology will automatically sync your brush head with the White+ mode to help whiten your teeth.

Built-in pressure sensor



You might not notice if you're brushing too hard, but your DiamondClean 9000 will. If you need to ease up, the toothbrush will make a pulsing sound. It's a heads up to let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

Long-life battery



Your DiamondClean 9000 will last for up to two weeks on a single charge with regular use.

DiamondClean 9000

Sonic electric toothbrush with app

HX9911/76

Specifications

Connectivity

Bluetooth® wireless technology: Connected brushing app

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): 14 days***

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: White

Service

Warranty: 2-year limited warranty

Compatibility

Android compatibility: Android phones, Bluetooth 4.0 enabled tablets

iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Replacement reminder: To always ensure best results, reminder icon lights up

Items included

Handle: 1 DiamondClean 9000

Travel case: 1

Brush heads: 1 G3 Premium Gum Care

Glass charger: 1

Charger base: 1

Cleaning performance

Health benefits: Up to 7x healthier gums*

Whitening benefits: Whitens teeth in 1 day****

Timer: BrushPacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday clean

3 intensities: Low, Medium, High

Deep Clean+: For an invigorating deep clean

Gum Health: Special attention to molars

White+: To remove surface stains

Smart sensor technology

BrushSync Replacement Reminder: Always know when to, replace brush heads

Pressure sensor: Alerts when brushing too hard

BrushSync Technology: Connects smart handle and, smart brush head

