

PHILIPS

sonicare

Sonic electric
toothbrush

CleanCare

1 mode

1 brush head

HX5350



Better plaque removal

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

Proven to improve oral health

- Clinically proven safe and gentle

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time

Provides a superior clean

- Sonicare dynamic sonic action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

Offers a customised brushing experience

- Easy-start programme to acclimate to the Sonicare experience

Highlights

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

Angled brush head neck



The uniquely angled neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Safe and gentle




Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

Smartimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience

Specifications

Technical specifications

Battery: Rechargeable
Charging time: 24 hour(s)
Power consumption: 0.9/1.1 W

Items included

Charger base: 1

Cleaning performance

Brush heads: 1 e Series standard brush head
Speed: Up to 31,000 brush movements/min
Timer: Smartimer

