

PHILIPS

AVENT

Manual breast pump set with steriliser

Pump, feed and sterilise

Includes 3-in-1 steriliser

Comfort Manual breast pump

Bottles, brush and soothers

SCD293/00



More comfort, more milk* Pump, feed and sterilise

Manual breast pump set with steriliser

The breastfeeding set is ideal for mums who want a full package, from pumping to feeding and cleaning up. It includes a variety of accessories: a manual breast pump, steriliser, soothers and bottles.

Feed

- Includes our Natural bottle and teat for natural latch-on

3-in-1 steriliser

- Adjustable size
- Natural steam sterilising kills 99.9% of harmful germs
- Suitable for wide neck and narrow neck bottles

Express - Manual Breast Pump

- Relaxing position thanks to its special design
- Soft massage cushion gently stimulates milk flow
- Compact lightweight design
- Few separate parts and intuitive design

Highlights

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a soft, velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let-down.

Easy on-the-go



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on the go more discreet.

Easy to set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

Includes our Natural bottle



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

3-in-1 steriliser



Due to its adjustable size, this 3-in-1 Philips Avent steriliser perfectly fits the items you would like to sterilise and always takes up the smallest amount of kitchen space. The 3-in-1 functionality gives you: 1) Small size configuration for sterilising soothers 2) Medium size configuration for sterilising breast pumps, toddler plates, knives and forks 3) Large size configuration for sterilising 6 baby

Natural steam sterilising

The steriliser uses natural steam to sterilise baby bottles and other products, killing 99.9% of harmful germs without using chemicals. Perfect for giving you peace of mind that all your baby bottles and other products are sterile.

Fits all bottle types



The steriliser has been designed to fit both wide neck and narrow neck bottles as well as different breast pumps and accessories. This gives you the freedom to choose the bottles and accessories that best suit you and your baby. At full capacity the steriliser holds six 11 oz/330 ml Philips Avent bottles.

Specifications

Power

Voltage: 220 V (3-in-1 electric steam steriliser)

Material

Bottle: Polypropylene, BPA-free*

Breast pump: BPA-free*, Polypropylene

What is included

Silicone Soother: 3 pcs

Breast pump body with handle: 1 pcs

Tongs: 1 pcs

Sealing disc for milk storage: 1 pcs

Extra soft Newborn Flow teat: 2 pcs

Natural bottle 125 ml/4 oz: 2 pcs

Snap on hygienic cap: 3 pcs

Natural bottle 9 oz: 2 pcs

Bottle and teat brush: 1 pcs

3-in-1 steriliser: 1 pcs

Breast pad sample packs: 1 pack (2 pcs)

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Development stages

Stage: 0–6 months

