

# IHE Integration Statement forAudit



Vendor	Product Name	Version	Date
Forcare	forAudit	2022-3	2022-08-16

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client (ITI-1)	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20) Audit Record Repository (ITI-20)	
Retrieve ATNA Audit Event (ATNA)	Audit Record Repository (ITI-81)	
Cross-Enterprise User Assertion (XUA)	X-Service-User (ITI-40)	Subject-Role
Internet User Authorization (IUA)	Authorization Client (ITI-71, ITI-72) Authorization Server (ITI-71) Resource Server (ITI-72)	JWT Token SAML Token JWT Token SAML Token
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer (ITI-18, ITI-43)	
Enterprise User Authentication (EUA)	Client Authentication Agent (ITI-2, ITI-3, ITI-4) Kerberized Server (ITI-4)	

Internet address for vendor's IHE information [www.forcare.com](http://www.forcare.com)

## Links to Standards Conformance Statements for the Implementation

HL7	<a href="http://www.forcare.com">www.forcare.com</a>
DICOM	<a href="http://www.forcare.com">www.forcare.com</a>

## Links to general information on IHE

In North America: [www.ihe.net](http://www.ihe.net)      In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)      In Japan: [www.ihe-j.org](http://www.ihe-j.org)

© 2006-2022 Forcare Holding B.V. All rights reserved.