

IHE Integration Statement forIndex



Vendor	Product Name	Version	Date
Forcare	forIndex	2022-3	2022-08-16

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client (ITI-1)	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
Cross-Enterprise User Assertion (XUA)	X-Service-Provider (ITI-40)	Subject-Role PurposeOfUse
	X-Service-User (ITI-40)	Subject-Role PurposeOfUse
Internet User Authorization (IUA)	Resource Server (ITI-72)	SAML Token
Cross-Enterprise Document Sharing (XDS.b)	Document Registry (ITI-18, ITI-42)	Patient Identity Feed (HL7 V2) (ITI-8)
		Document Metadata Update (ITI-57)
		On-Demand Documents (ITI-61) Reference ID
Remove Metadata and Documents (RMD)	Document Registry (ITI-18, ITI-62)	
Multi-Patient Queries (MPQ)	Document Registry (ITI-51)	
Patient Demographics Query (PDQ)	Patient Demographics Supplier (ITI-21)	
Patient Demographics Query HL7 V3 (PDQV3)	Patient Demographics Supplier (ITI-47)	Continuation Option
Patient Identifier Cross-referencing for MPI (PIX)	Patient Identifier Cross-reference Manager (ITI-8, ITI-9)	
Patient Identifier Cross-referencing HL7 V3 (PIXV3)	Patient Identifier Cross-reference Manager (ITI-45)	
Document Metadata Subscription (DSUB)	Document Metadata Notification Broker (ITI-52, ITI-53)	Document Metadata Publish Recipient (ITI-54)
	Document Metadata Publisher (ITI-54)	

Internet address for vendor's IHE information www.forcare.com

Links to Standards Conformance Statements for the Implementation

HL7	www.forcare.com
DICOM	www.forcare.com

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--

© 2006-2022 Forcare Holding B.V. All rights reserved.