

# IHE Integration Statement forConnect



Vendor	Product Name	Version	Date
Forcare	forConnect	2022-3	2022-08-16
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below			
Integration Profile	Actor	Option	
Consistent Time (CT)	Time Client (ITI-1)		
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)		
Cross-Enterprise User Assertion (XUA)	X-Service-User (ITI-40)		
Internet User Authorization (IUA)	Resource Server (ITI-72)	SAML Token	
Cross-Enterprise Document Sharing (XDS.b)	Document Source (ITI-41)	Document Addendum	
		Document Replacement	
	Document Consumer (ITI-18, ITI-43)	Document Transformation	
		Document Metadata Update	
Cross-Enterprise Document Reliable Interchange (XDR)	Document Source (ITI-41)	On-Demand Documents	
		Delayed Document Assembly	
Cross-Enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Source (ITI-15)	Update Metadata (ITI-57)	
		Transmit Home Community ID	
Cross-Enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Source (ITI-15)	Set of DICOM Instances	
		PDF Report	
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Source (ITI-41, RAD-68)	Text Report	
		Imaging Document Consumer (RAD-55)	
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Source (ITI-41, RAD-68)	Set of DICOM Instances (RAD-55, RAD-69)	
		PDF Report	
	Imaging Document Consumer (RAD-55, RAD-69)	CDA Wrapped Text Report	
		CDA Imaging Report with Structured Headings	
Cross Community Access (XCA)	Responding Gateway (ITI-38, ITI-39)		
Patient Identifier Cross-referencing for MPI (PIX)	Patient Identifier Cross-Reference Consumer (ITI-9)		
Patient Identifier Cross-referencing HL7V3 (PIXV3)	Patient Identifier Cross-Reference Consumer (ITI-45)		
Basic Patient Privacy Consents (BPPC)	Content Creator		
Scanned Documents (XDS-SD)	Content Creator		
Cross-Enterprise Document Workflow Profile (XDW)	Content Creator	XDW Cross-Community with XDR	
	Content Consumer	Document Import Option	
	Content Updater	Document Import Option	
Cross Enterprise Basic eReferral Workflow Definition Profile (XBeR-WD)	Referral Requester	XDW Cross-Community with XDR	
Internet address for vendor's IHE information <a href="http://www.forcare.com">www.forcare.com</a>			

## Links to Standards Conformance Statements for the Implementation

HL7

[www.forcare.com](http://www.forcare.com)

---

DICOM

[www.forcare.com](http://www.forcare.com)

---

**Links to general information on IHE**

---

In North America: [www.ihe.net](http://www.ihe.net)

In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)

In Japan: [www.ihe-j.org](http://www.ihe-j.org)

---

© 2006-2022 Forcare Holding B.V. All rights reserved.

---